

Importance of Acrylamide Determination


Acrylamide is an essential industrial chemical with an estimated worldwide production of 200 million Kg/year. It is used as a grouting agent in construction, a paper making aid, a soil conditioning agent, in ore processing, in sewage treatment, and as an additive (coagulant) in water treatment. Also a component of cigarette smoke, Acrylamide is a known carcinogen in laboratory animals, impairing fertility in male animals and causing nerve damage in humans with industrial exposure.

In addition to its industrial uses, Acrylamide is also found as a natural byproduct of the cooking process. Methods in which temperature exceeds 120 °C, such as baking, frying, grilling, and toasting can cause the amino acid asparagine (found in certain foods) to react with reducing sugars such as glucose (via the Maillard reaction) to produce Acrylamide. High carbohydrate foods that are baked or fried at high temperatures contain the highest levels of Acrylamide. Acrylamide is not found in raw or boiled foods that are high in carbohydrates or in meats.

The Abraxis Acrylamide ELISA allows the determination of up to 40 samples in duplicate determination. Only a few grams or milliliters of sample are required. The test requires an extraction and clean-up step, and derivatization prior to assaying. The actual ELISA can be performed in less than 1.5 hours.

Performance Data

Test sensitivity: The limit of quantitation for Acrylamide calculated as $X_n \pm 3SD$ ($n=20$) in the various matrices is as follows: Aqueous solutions 2.5 ppb (lower LOD can be obtained with SPE); Food samples 25 ppb (contact Abraxis for alternative procedures to obtain LOD <25 ppb in food samples).



Determinations closer to the middle of the calibration range of the test yield the most accurate results.

Test reproducibility: Coefficients of variation (CVs) for standards: <10%; CVs for samples: <20%.

Selectivity: This ELISA recognizes Acrylamide and not related compounds:

Cross-reactivities:	Acrylamide	100%
	Acrylonitrile	2.1%

The following compounds demonstrated no reactivity in the Abraxis Acrylamide ELISA at concentrations up to 10,000 ppb: Acrolein, Acrylic acid, Asparagine, Aspartamine, Aspartic acid, Glutamic acid, Glutamine, Methacrylamide, Methyl acrylate, 2-Pyrrolidone, 2-Pyrroglutamic acid.

Samples: To eliminate matrix effects from food samples, a sample clean-up is required. See Preparation of Samples (Section H).

General Limited Warranty: Abraxis LLC warrants the products manufactured by the Company against defects and workmanship when used in accordance with the applicable instructions for a period not to extend beyond the product's printed expiration date. **Abraxis makes no other warranty, expressed or implied. There is no warranty of merchantability or fitness for a particular purpose.**

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, www.atzlabs.com; www.lifetechindia.com

Acrylamide-ES ELISA, Microtiter Plate

Enzyme-Linked Immunosorbent Assay for the Determination

of Acrylamide in Food Samples

Product No. 515680

1. General Description

The Acrylamide ELISA is an immunoassay for the detection of Acrylamide. This test is suitable for the quantitative and/or qualitative detection of Acrylamide in food samples. Samples requiring action should be confirmed by HPLC, GC/MS, or other conventional methods.

2. Safety Instructions

The standard solutions in this test kit contain small amounts of Acrylamide. In addition, the substrate solution contains tetramethylbenzidine and the stop solution contains diluted sulfuric acid. Avoid contact of stopping solution with skin and mucous membranes. If these reagents come in contact with the skin, wash with water.

3. Storage and Stability

The Acrylamide ELISA Kit should be stored in the refrigerator (4–8°C). The plate, standard/control, color and stop solutions must be allowed to reach room temperature (20–25°C) before use. Reagents may be used until the expiration date on the box. The Enzyme Conjugate is supplied in a lyophilized form (3 vials). Before each assay, the required volume of lyophilized conjugate must be reconstituted with the appropriate diluent (see Test Preparation section). Reconstitute only the amount needed for the samples to be run, as the reconstituted solution will only remain viable for one day. Standards, controls and samples need to be derivatized before assaying and **must be derivatized and run the same day.**

4. Test Principle

The test is a direct competitive ELISA based on the recognition of Acrylamide by specific antibodies. Acrylamide, when present in a sample, and an Acrylamide-enzyme conjugate compete for the binding sites of rabbit anti-Acrylamide antibodies in solution. The Acrylamide antibodies are then bound by a second antibody (goat anti-rabbit) immobilized on the microtiter plate. After a washing step and addition of the substrate solution, a color signal is produced. The intensity of the blue color is inversely proportional to the concentration of Acrylamide present in the sample. The color reaction is stopped after a specified time and the color is evaluated using an ELISA reader. The concentrations of the samples are determined by interpolation using the standard curve constructed with each run.

5. Limitations of the Acrylamide ELISA, Possible Test Interference

Numerous organic and inorganic compounds commonly found in samples have been tested and found not to interfere with this test. However, due to the high variability of compounds that might be found in samples, test interferences caused by matrix effects can not be completely excluded. Mistakes in handling the test can also cause errors. Possible sources for such errors can be:

Inadequate storage conditions of the test kit, incorrect pipetting sequence or inaccurate volumes of the reagents, too long or too short incubation times during the immune and/or substrate reaction, extreme temperatures during the test performance (lower than 2°C or higher than 8°C during the first incubation or lower than 10 °C or higher than 30°C during the color incubation step), inappropriate sample clean up or derivatization.

The Abraxis Acrylamide ELISA kit provides screening results. As with any analytical technique (GC, HPLC, etc.), samples requiring some action should be confirmed by an alternative method.

Working Instructions

A. Materials Provided

1. Microtiter plate coated with a second antibody (12 X 8 strips).
2. Acrylamide Standards/Calibrators (7): 0; 2.5; 5; 10; 25; 50, and 200 ng/mL, 3.0 mL each. Refer to Section I for derivatization procedure.
3. Acrylamide Control, 1 vial, 3.0 mL. Refer to Section I for derivatization procedure.
4. Acrylamide-HRP Enzyme Conjugate, 3 vials (lyophilized), 2 mL/vial after reconstitution.
5. Enzyme Conjugate Diluent, 8 mL.
6. Antibody Solution, rabbit anti-Acrylamide, 6 mL.
7. Derivatization Kit (PN 515676), containing: a) Assay Buffer, 100 mL. Use to neutralize samples after derivatization; b) Derivatization Reagent 2 X 2 mL each (lyophilized, see reagent preparation section); c) Sample Diluent Solution (60% Methanol in water), 2 X 20 mL, use to dilute samples with high concentrations (outside the calibration range of the assay); d) Derivatization Reagent reconstitution solution, 6 mL.
8. Wash Solution (5X) Concentrate, 100 mL.
9. Color (Substrate) Solution (TMB), 16 mL.
10. Stop Solution, 12 mL (contains diluted H₂SO₄; should be handled with care).

B. Test Preparation

Micro-pipetting equipment and pipette tips for pipetting the standards and the samples are necessary. We recommend using a multi-channel pipette or a stepping pipette for adding the conjugate, antibody, substrate and stop solutions in order to equalize the incubation periods of the solutions on the entire microtiter plate. Please use only the reagents and standards from one package lot in one test, as they have been adjusted in combination.

- Adjust the microtiter plate, standards/control, color, and stop solutions to room temperature before use. Antibody solution, conjugate diluent and HRP conjugate should be removed from refrigerator approximately 15 minutes prior to assay. **Wash solution needs to be cold.**
- Remove the number of microtiter plate strips required from the foil bag. The remaining strips are stored in the foil bag and zip-locked closed. Store the remaining kit in the refrigerator (4-8°C).
- The standard solutions, antibody, substrate and stop solutions are ready to use and do not require any further dilutions.
- The Derivatization solution is lyophilized (2 vials). Before each assay, calculate the volume needed (when reconstituted, each vial will provide enough derivatization solution for approx. 80 wells). Reconstitute only the amount necessary for the samples to be analyzed. Once reconstituted, the derivatization solution will only remain viable for 1 week. If the assay requires >80 wells, a second vial will need to be prepared and combined with the first vial before use. To reconstitute, add 2 mL of Derivatization Diluent to each vial of Derivatization Solution, allow to sit for 5 minutes and vortex for 5-10 seconds.
- The enzyme conjugate provided is lyophilized (3 vials). Before each assay, calculate the volume of conjugate needed (when reconstituted, each vial will provide enough conjugate for approx. 40 wells). Reconstitute only the amount necessary for the samples to be analyzed. Once reconstituted, the conjugate solution will only remain viable for 1 day. If the assay requires >40 wells, a second vial of conjugate will need to be prepared and combined with the first vial before use. To reconstitute, add 2 mL of Enzyme Conjugate Diluent to each vial of Enzyme Conjugate, allow to sit for 5 minutes and vortex for 5-10 seconds, store in dark prior to assay.
- Dilute the Wash Buffer concentrate at a ratio of 1:5. If using the entire bottle (100 mL), add to 400 mL of deionized or distilled water. Store at 4-8°C, **remove from refrigerator just prior to washing plates.**

C. Assay Procedure

- Add 50 µL of the derivatized standard solutions, control, and derivatized samples or derivatized sample extracts (Please see section H. Preparation of samples and section I. Derivatization Procedure) into the wells of the test strips according to the working scheme given. We recommend using duplicates or triplicates.
- Add 50 µL of the reconstituted **enzyme conjugate** solution to the individual wells successively using a multi-channel pipette or a stepping pipette.
- Add 50 µL of **antibody solution** to the individual wells successively using a multi-channel pipette or a stepping pipette. Cover the wells with parafilm or tape and mix the contents by moving the strip holder in a circular motion on the benchtop for about 60 seconds. Be careful not to spill contents.
- Incubate the strips for 60 minutes at **refrigerated temperature** (2°C to 8°C).
- After incubation, remove the covering and vigorously shake the contents of these wells into a sink. Wash the strips **three times** using the **cold 1X** washing buffer solution. Use at least a volume of 250 µL of washing buffer for each well and each washing step. Remaining buffer in the wells should be removed by patting the plate dry on a stack of paper towels.
- Add 150 µL of **substrate (color) solution** to the wells. Cover the wells with parafilm or tape and mix the contents by moving the strip holder in a circular motion on the benchtop for **60-90 seconds**. Incubate the strips for 20-30 minutes at **room temperature**. Protect the strips from direct sunlight.
- Add 100 µL of **stop solution** to the wells in the same sequence as for the substrate solution and mix the contents by moving the strip holder in a circular motion on the benchtop for 15-30 seconds.
- Read the absorbance at 450 nm using a microplate ELISA photometer within 15 minutes after the addition of the stopping solution.

D. Evaluation

The evaluation of the ELISA can be performed using commercial ELISA evaluation programs (4-Parameter (preferred) or Logit/Log). For manual evaluation, calculate the mean absorbance value for each of the standards. Calculate the %B/B₀ for each standard by dividing the mean absorbance value for each standard by the Zero Standard (Standard 0) mean absorbance. Construct a standard curve by plotting the %B/B₀ for each standard on the vertical linear (y) axis versus the corresponding Acrylamide concentration on the horizontal logarithmic (x) axis on graph paper. %B/B₀ for samples will then yield levels in ppb of Acrylamide by interpolation using the standard curve after multiplying by a factor of 10. Samples showing lower concentrations of Acrylamide compared to Standard 1 (2.5 ng/mL) are considered as negative. Samples showing a higher concentration than Standard 6 (200 ng/mL) must be diluted further to obtain accurate results.

E. Additional Materials (not included with the test kit)

- | | |
|---|---|
| 1. Micro-pipettes with disposable plastic tips (10-200 and 200-1000 µL) | 9. Dual Incubator (2-8 °C and 50 °C) or equivalent |
| 2. Multi-channel pipette (10-250 µL) or stepper pipette with plastic tips (10-250 µL) | 10. Timer |
| 3. Microtiter plate reader (wave length 450 nm) | 11. Tape or Parafilm |
| 4. Microcentrifuge, capable of spinning up to 15,000 X g | 12. Serological pipettes, 20 mL capacity or repipet |
| 5. Vortexer | 13. SPE array manifold with vacuums source |
| 6. Orbital Mixer | 14. 2.5 mL collection plate (96 wells) |
| 7. 50 mL sample extraction tube with filter; 4 mL glass vials with caps | 15. 2 mL microcentrifuge tube capable of spinning up to 15,000 X g |
| 8. Methanol | 16. SPE columns (2): Tube # 1 array column, 2 mL, cat # 395010. Tube # 2 array column, 2 mL, cat # 395009 |

F. Additional Reagents (not included with the test kit)

Methanol, 60% Methanol in water, distilled or deionized Water.

G. Working Scheme

The microtiter plate consists of 12 strips of 8 wells, which can be used individually for the test. The standards must be run with each test. Never use the values of standards which have been determined in a test performed previously.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 4	Samp 1									
B	Std 0	Std 4	Samp 1									
C	Std 1	Std 5	Samp 2									
D	Std 1	Std 5	Samp 2									
E	Std 2	Std 6	etc.									
F	Std 2	Std 6	etc.									
G	Std 3	Control										
H	Std 3	Control										

Std 0-Std 6: Standards

0; 2.5; 5; 10; 25; 50, 200 ppb

Control: 20 +/- 5 ppb

Samp1, Samp2, etc.: Samples

H. Preparation of Samples

Potato Chip/Corn Chip Sample Extraction/Clean Up (Alternate clean-up procedures are available)

- Crush sample using mortar and pestle.
- Weigh 2.0 g of a representative sample into a 50 mL sample extraction tube.
- Add 20 mL of distilled or deionized water.
- Place sample in an orbital mixer and mix for 30 minutes.
- Remove tube from rotary mixer and place on table top for approx. 5 minutes to allow sedimentation.
- Using the filter plunger, push filter into sample extract to expedite separation.
- Transfer 4 mL of the filtered sample extract into 2 separate 2 mL microcentrifuge tubes (2 mL each)
- Centrifuge tube for 5 minutes at 10,000-15,000 X g at room temperature.
- Carefully transfer 1.0 mL of supernatant from each of the 2 tubes into one clean sample vial for further clean up (be sure to avoid removing top oily layer and bottom solid layer). This combined sample is used in step 11.
- Condition a Multi Mode array column with 1.5 mL methanol, followed by 4.0 mL of deionized water. Both eluted conditioning solutions are discarded.
- A 1.8 mL aliquot of the extract from step 9 is passed through the column at a flow rate of 0.5-1.0 mL/minute. Eluent is collected into a clean 2.5 mL collection plate.
- Rinse the column with 0.6 mL of deionized water at a flow rate of 0.5-1.0 mL/minute. Collect the rinse in the same well of the collection plate containing the previously eluted extract (step 11).
- Condition an ENV+ array column with 2.0 mL of Methanol, followed by 2.0 mL of deionized water. Both eluted conditioning solutions are discarded.
- Pass the eluant (Step H. 11) and rinse (Step H. 12) through the column at a flow rate of 0.5-1.0 mL/minute. Discard eluent.
- Rinse the column with 2 X 1.0 mL of deionized water at a flow rate of 0.5-1.0 mL/minute. Discard rinse.
- Allow column to dry completely (approx. 15 min) and then elute acrylamide from the column with 2 X 0.9 mL of 60% Methanol/deionized water at a flow rate of 0.5-1.0 mL/minute. Collect eluent into a clean 2.5 mL collection plate. This eluent is derivatized in steps I. 1-5.

I. Derivatization Procedure

- 250 µL of sample extract (Step H.16) is added to a labeled glass vial with screw cap.
- Add 50 µL of derivatization reagent to each sample.
- Vortex vigorously for 10-15 seconds.
- Incubate at 47-53 °C for 60 minutes. Allow sample to cool for 15 minutes.
- Add 2.0 mL of Acrylamide Assay Buffer to cooled derivatized sample
- Analyze as sample (Assay Procedure, step 1). Derivatization and analysis must be performed the same day.

The ELISA result must be multiplied by a factor of 10 to obtain the final Acrylamide concentration in the sample (the multiplication factor is necessary to account for the sample extraction/dilution). Highly contaminated samples (those outside of the calibration range of the assay) must be diluted and re-analyzed.

For additional extraction procedures (other matrices), and lower LOD (<25 ppb) please contact Abraxis LLC.