

ALK Break Apart FISH Probe Kit

Introduction

The ALK Break Apart FISH Probe Kit is designed to detect rearrangements in the human *ALK* gene located on chromosome band 2p23.2. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other ALK aberrations such as deletions or amplifications. Initially discovered in anaplastic large cell lymphoma (ALCL), rearrangements of *ALK* – also known as *CD246* or *NBLST3* – have since been found in many types of malignancies, including B- and T-cell lymphomas, plasmacytomas, neuroblastoma, esophageal, breast, kidney, colon thyroid, lung and other cancers. A significant percentage of non-small cell lung cancer (NSCLC) cases harbor *ALK* gene abnormalities.

Intended Use

To detect rearrangements in the human *ALK* gene located on chromosome band 2p23.2.


Cont.

LSP ALK 5' FISH Probe
LSP ALK 3' FISH Probe

Color

CytoGreen
CytoOrange

Probe Design



Not to Scale

LSP ALK 5' FISH Probe covers the 5' and the center sequences of the *ALK* gene. LSP ALK 3' FISH Probe covers the 3' end and the neighboring downstream region. The two probes are designed to recognize sequences on both sides of a common breakpoint that is located inside the *ALK* gene.

Cat. No.

CT-PAC009-10-GO

Volume

10 Tests (100 µL)

Signal Pattern Interpretation

Normal Pattern


2OG*

Abnormal Pattern

Other Patterns

*Overlapping orange and green signals can appear as yellow.

- 1) Minoo P & Wang HY. *Int J Clin Exp Pathol.* 5(5):397-410 (2012).
- 2) Chiarle R, et al. *Nat Rev Cancer.* 8(1):11-23 (2008).
- 3) Salido M, et al. *J Thorac Oncol.* 6(1):21-7 (2011).
- 4) Kwak EL, et al. *N Engl J Med.* 363(18):1693-703 (2010).
- 5) Sasaki T, et al. *Eur J Cancer.* 46(10):1773-80 (2010).



CytoTest Inc.
9430 Key West Ave., Suite 210
Rockville, MD 20850, USA

* CE IVD only available in certain countries. All other countries are either ASR or RUO. Please contact your local dealer or our headquarters for more information.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com