

To order contact:  
306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222  
Mobile: +91-9810521400 Fax: +91-11-42208444  
Email: [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)  
Web: [www.atzlabs.com](http://www.atzlabs.com)



## Sample Stabilization and Recovery Quick Reference Protocol 1.5-ml Screw-Cap Tube and 96-well Plate

RNAstable™ protects RNA samples from degradation at room temperature. Each tube or plate contains RNAstable as a coating at the bottom of the tube or well that protects picogram to microgram amounts of RNA. This medium is completely dissolvable and ensures total sample recovery.

### Stabilize for Storage

- Open cap on tube or remove seal from 96-well plate.
- Add up to 100  $\mu\text{l}$  [ $\leq 100 \mu\text{g}$ ] of the RNA sample directly into each tube or well.
- For convenient air-drying of volumes  $\leq 20 \mu\text{l}$ , leave tube open overnight in laminar flow hood.
- For volumes  $\geq 20 \mu\text{l}$ , or for accelerated drying of all volumes, please SpeedVac without heat.
- Once dried, cap tube or seal 96-well plate. Store at room temperature (15-25°C) and protect from moisture. Either:
  - 1) Store in a dry storage cabinet.
  - 2) Heat seal the moisture-barrier bag, enclosing desiccant pack.
- Dry storage cabinets and additional moisture-barrier bags are available at [www.biomatrica.com](http://www.biomatrica.com).

### SpeedVac® Drying Time

Sample Volume	Drying Time
10-20 $\mu\text{l}$	30 minutes
20-30 $\mu\text{l}$	1 hour
30-100 $\mu\text{l}$	1.5 hours

For sample recovery, see reverse.

For more information, please refer to the RNAstable handbook at [www.lifetechindia.com/RNA-Stable-Detailed-Handbook.pdf](http://www.lifetechindia.com/RNA-Stable-Detailed-Handbook.pdf)

To order contact:  
306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222  
Mobile: +91-9810521400 Fax: +91-11-42208444  
Email: [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)  
Web: [www.atzlabs.com](http://www.atzlabs.com)



Sample Stabilization and Recovery  
Quick Reference Protocol  
1.5-ml Screw-Cap Tube  
and 96-well Plate

## Sample Recovery: ***Just Add Water***

- Add 10-100  $\mu$ l of H<sub>2</sub>O or other liquid to the tube or well containing stored sample.
- Hydrate for 15 minutes.
- Pipette gently to ensure complete mixing.

## Samples can be used directly in downstream applications:

- Quantitative Real-Time PCR
- Bioanalyzer and microarray analysis
- End-point PCR and gel analysis
- cDNA synthesis
- Reverse transcription

For more information, please refer to the RNastable handbook at [www.lifetechindia.com/RNA-Stable-Detailed-Handbook.pdf](http://www.lifetechindia.com/RNA-Stable-Detailed-Handbook.pdf)

**To order contact:**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400 Fax: +91-11-42208444

Email: [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)

Web: [www.atzlabs.com](http://www.atzlabs.com)