

(G418, Hygromycin B – HygroGold, Puromycin, Zeocin, Blasticidin S, Phleomycin)

The chemical structure shows a complex oligosaccharide composed of several sugar units. It includes a hexose unit with a methyl group and a hydroxyl group, linked to another hexose unit which has an amino group. This is further linked to a pyranose unit with an amino group, which is then linked to a furanose unit with a hydroxyl group. The structure is highly branched and contains multiple stereocenters indicated by wedges and dashes.


Selection in mammalian cells is usually achieved in three to seven days with concentrations ranging from 400 to 1000  $\mu\text{g}/\text{ml}$ . Cells that are dividing are affected sooner than those that are not.

CAS Number: 108321-42-2

- [pSELECT-neo-mcs](#), [pSELECT-neo-LacZ](#) : pSELECT plasmids offer all the features necessary to express a gene of interest at high levels in a large number of cell types.
- [pMONO-neo-mcs](#), [pMONO-neo-qfp](#) : pMONO plasmids are specifically designed for strong and constitutive expression of a gene of interest in a wide variety of cell lines.
- [pVITRO1-neo-mcs](#), [pVITRO1-neo-GFP/LacZ](#) : pVITRO allow the ubiquitous and constitutive co-expression of two genes of interest.
- [pVITRO2-neo-mcs](#), [pVITRO2-neo-GFP/LacZ](#) : pVIVO allow strong and sustained co-expression of two genes of interest in many tissues and organs.
- [psiRNA-h7SK G1 Neo](#) : psiRNA-h7SK offer all the features necessary to express shRNA under the human 7SK RNA polymerase III promoter.
- [pCpGfree-vitroNmcs](#), [pCpGfree-vltroNLacZ](#) : pCpGfree-vitro is a family of expression vectors completely devoid of CpG dinucleotides that are selectable in mammalian cells

G418 is shipped at room temperature and should be stored at -20°C. G418 is stable for at least 1 year when stored at -20°C.

## Hygromycin B - HygroGold™



Hygromycin B is an aminoglycoside antibiotic produced by *Streptomyces hygroscopicus*. **It inhibits protein synthesis** by interfering with translocation and causing mistranslation at the 70S ribosome. Hygromycin B is effective **on most bacteria, fungi and higher eukaryotes**.

Resistance to hygromycin is conferred by the **hph gene** from *E. coli*.

Hygromycin B is normally used at a concentration of 50-200 µg/ml in mammalian cells and 100 µg/ml in bacteria.

Two grades of Hygromycin B are available:

**Hygromycin B** (purity >85%)

**HygroGold™** (purity >98%)

Formula: C<sub>20</sub>H<sub>37</sub>N<sub>3</sub>O<sub>13</sub>, HCl

Molecular weight: 527.52

CAS Number : 31282-04-9

**Bulk quantities available at highly discounted prices :** Contact [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)

## **Cloning vectors using Hygromycin B**

InvivoGen offers several selection vectors using Hygromycin as selection antibiotic:

- [pSELECT-hygro-mcs](#), [pSELECT-hygro-LacZ](#) : pSELECT plasmids offer all the features necessary to express a gene of interest at high levels in a large number of cell types.
- [pMONO-hygro-mcs](#), [pMONO-hygro-gfp](#) : pMONO plasmids are specifically designed for strong and constitutive expression of a gene of interest in a wide variety of cell lines.
- [pVITRO1-hygro-mcs](#), [pVITRO1-hygro-GFP/LacZ](#) : pVITRO allow the ubiquitous and constitutive co-expression of two genes of interest.
- [pVITRO2-hygro-mcs](#), [pVITRO2-hygro-GFP/LacZ](#) : pVIVO allow strong and sustained co-expression of two genes of interest in many tissues and organs.
- [psiRNA-h7SK G1 Hygro](#) : psiRNA-h7SK offer all the features necessary to express shRNA under the human 7SK RNA polymerase III promoter.
- [pCpGfree-vitroHmcs](#), [pCpGfree-vitroHLacZ](#) : pCpGfree-vitro is a family of expression vectors completely devoid of CpG dinucleotides that are selectable in mammalian cells

## **Contents and Storage**

Hygromycin B and HygroGold™ are provided as 100 mg/ml yellow solutions. HygroGold™ is also provided as a powder. Products are shipped at room temperature. Store at -20°C. Hygromycin B solutions are stable for at least one year when stored at -20°C.



## Puromycin

Puromycin is an aminonucleoside antibiotic produced by *Streptomyces alboniger*. It specifically **inhibits peptidyl transfer** on both prokaryotic and eukaryotic ribosomes. This antibiotic inhibits the growth of **Gram positive bacteria** and **various animal and insect cells**. Puromycin can also be used in some particular conditions for the selection of *E. coli* transformants.

Resistance to puromycin is conferred by the **Pac gene** encoding a puromycin N-acetyl-transferase (PAC) that was found in a *Streptomyces* producer strain.

Animal cells are generally sensitive to concentrations from 1 to 10 µg/ml.

Formula: C<sub>22</sub>H<sub>29</sub>N<sub>7</sub>O<sub>5</sub>, 2HCl

Molecular weight: 471.51

CAS Number: 58-58-2

**Bulk quantities available at highly discounted prices :** Contact [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)

## **Cloning vectors using Puromycin**

InvivoGen offers selection vectors using Puromycin as selection antibiotic:

- [pSELECT-puro-mcs](#), [pSELECT-puro-LacZ](#) : pSELECT plasmids offer all the features necessary to express a gene of interest at high levels in a large number of cell types

## **Contents and Storage**

Puromycin hydrochloride is provided as a colorless solution at 10 mg/ml.

Puromycin is shipped at room temperature and should be stored at -20°C. Puromycin is stable for one year when stored at -20°C.



## Zeocin™

InvivoGen is the sole worldwide producer of Zeocin™, a copper-chelated glycopeptide antibiotic produced by *Streptomyces CL990*.

Zeocin™ causes cell death by intercalating into DNA and cleaving it. This antibiotic **is effective on most aerobic cells** and is therefore useful for selection in **bacteria, eukaryotic microorganisms, plant and animal cells**. Resistance to Zeocin™ is conferred by the **Sh ble** gene product which inactivates Zeocin™ by binding to the antibiotic.

Zeocin™ is used at a concentration of 50-300 µg/ml for selection in mammalian cells and 25 µg/ml for bacterial selection.

Empirical formula: C<sub>55</sub>H<sub>85</sub>O<sub>21</sub>N<sub>20</sub>S<sub>2</sub>Cu - HCl

Molecular weight: 1525

**Bulk quantities available at highly discounted prices :** Contact [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)

Cloning vectors using Zeocin™


Invivogen offers several selection vectors using Zeocin™ as selection antibiotic:

- [pSELECT-zeo-mcs](#), [pSELECT-gfpzeo-mcs](#), [pSELECT-zeo-LacZ](#), [pSELECT-gfpzeo-LacZ](#) : pSELECT plasmids offer all the features necessary to express a gene of interest at high levels in a large number of cell types.
- [pMONO-zeo-mcs](#), [pMONO-zeo-gfp](#) : pMONO plasmids are specifically designed for strong and constitutive expression of a gene of interest in a wide variety of cell lines.
- [psiRNA-h7SK G1 Zeo](#), [psiRNA-h7SK G1 GFPzeo](#), [psiRNA-DUO](#) : psiRNA-h7SK offer all the features necessary to express shRNA under the human 7SK RNA polymerase III promoter.
- [pCpG-free](#): pCpGfree is an expression vectors completely devoid of CpG dinucleotides that is selectable in mammalian cells

Contents and Storage

Zeocin™ is provided as a blue solution at 100 mg/ml.

## **Blasticidin S**



Blasticidin S is a peptidyl nucleoside antibiotic isolated from *Streptomyces griseochromogenes*. It specifically **inhibits protein synthesis in both prokaryotes and eukaryotes** by interfering with the peptide bond formation in the ribosomal machinery.

Resistance to blasticidin is conferred by the blasticidin resistance gene from *Bacillus cereus* (**bsr**) which codes for blasticidin deaminase.

Typically, bacteria are sensitive to blasticidin concentrations of 25-100 µg/ml, and mammalian cells to 1-10 µg/ml.

Formula: C<sub>17</sub>H<sub>26</sub>N<sub>8</sub>O<sub>5</sub>, 1 HCl

Molecular weight: 458.5

CAS Number: 2079-00-7

**Bulk quantities available at highly discounted prices :** Contact [customerservice@atclabs.com](mailto:customerservice@atclabs.com)

Cloning vectors using Blasticidin S

InvivoGen offers several selection vectors using Blasticidin S as selection antibiotic:

- [pSELECT-blasti-mcs](#), [pSELECT-blasti-LacZ](#) : pSELECT plasmids offer all the features necessary to express a gene of interest at high levels in a large number of cell types.
- [pMONO-blasti-mcs](#), [pMONO-blasti-gfp](#) : pMONO plasmids are specifically designed for strong and constitutive expression of a gene of interest in a wide variety of cell lines.
- [pVITRO1-blasti-mcs](#), [pVITRO1-blasti-GFP/LacZ](#) : pVITRO allow the ubiquitous and constitutive co-expression of two genes of interest.
- [pVITRO2-blasti-mcs](#), [pVITRO2-blasti-GFP/LacZ](#) : pVITRO allow strong and sustained co-expression of two genes of interest in many tissues and organs.
- [psiRNA-h7SK G1 Blasti](#) : psiRNA-h7SK offer all the features necessary to express shRNA under the human 7SK RNA polymerase III promoter.
- [pCpGfree-vitroBLacZ](#) : pCpGfree-vitro is a family of expression vectors completely devoid of CpG dinucleotides that are selectable in mammalian cells

Contents

Blasticidin is provided as a colorless solution at 10 mg/ml. Blasticidin is shipped at room temperature and should be stored at -20°C. Blasticidin is stable up to 1 year when stored at -20°C.



## **Phleomycin**

Phleomycin is a glycopeptide antibiotic of the bleomycin family, isolated from a mutant strain of *Streptomyces verticillus*. It binds and intercalates DNA thus destroying the integrity of the double helix. Phleomycin is **active against most bacteria, filamentous fungi, yeast, plant and animal cells.**

Use of phleomycin is recommended for cells poorly sensitive to Zeocin™, i.e. filamentous fungi and some yeasts. Phleomycin resistance is conferred by the **Sh ble** gene from *Streptoalloteichus hindustanus* which encodes a protein that binds to phleomycin, inhibiting its DNA cleavage activity.

Typically, phleomycin is used at a concentration of 10 µg/ml for yeasts and 25-150 µg/ml for filamentous fungi.

Empirical formula: C<sub>55</sub>H<sub>85</sub>O<sub>21</sub>N<sub>20</sub>S<sub>2</sub>Cu-HCL

Molecular weight: 1525

CAS Number: 11006-33-0

**Bulk quantities available at highly discounted prices :** Contact [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)

## Contents and Storage

Phleomycin is provided as a blue solution at 20 mg/ml or as a powder.

Phleomycin is shipped at room temperature. Store the solution at -20°C and the powder at 4°C. Phleomycin is stable for at least 1 year when properly stored.