

GenePurgeDirect® DNA/RNA Releasing Agent

General Protocol

Version: 1.0
Revision date: 31-07-2014

Description

GenePurgeDirect® is composed of proprietary polymeric materials that quickly facilitate the release of nucleic acids from cells in a form suitable for PCR. By segregating inhibitors that are released during lysis and preservation agents that may interfere with amplification, GenePurgeDirect® provides amplifiable nucleic acids from minute amounts of material.

Protocols

GenePurgeDirect® General Protocol

1. Place 1-5µl of whole blood / 1 isolated colony / 1ul of cells at 10^3 - 10^8 cells/ml into the bottom of a standard amplification for each specimen.
2. Resuspend the contents of the GenePurgeDirect® tube by inverting 10- 20 times or briefly vortexing.
3. Add 20µl of the resuspended GenePurgeDirect® to each tube, pulse vortex briefly to mix.
4. Perform the thermal cycler protocol (A) or microwave protocol (B).
5. Once program is completed, centrifuge sample tubes at 5,000 x g for 1 minute.
6. Transfer supernatant into a new tube for use as template for PCR. Use 1-10µl of supernatant per 20-100µl amplification reaction.
 - PCR can be performed directly in the GenePurgeDirect® treatment tube; add amplification reagents for a final volume of 100µl with all components at a final concentration of 1X.
7. Perform amplification reaction according to your optimized protocol.
 - The first cycle should have a denaturing time of 2-4 minutes at 94-95°C.
 - For some primer sets a re-titration of the magnesium concentration may be required.

(A) Thermal Cycler Lysis Protocol

Place samples onto thermal cycler, with a heated lid, with the following GenePurgeDirect® program:

Step	Temperature	Time
1.	65°C	30 sec.
2.	8°C	30 sec.
3.	65°C	90 sec.
4.	97°C	180 sec.
5.	8°C	60 sec.
6.	65°C	180 sec.
7.	97°C	60 sec.
8.	65°C	60 sec.
9.	80°C	hold

(B) Microwave Lysis Protocol

We have found that the microwave treatment of specimens affords a rapid sample preparation and facilitates the amplification of the more intractable types of specimens.

I. Evaluation of microwave

Perform the following experiment to determine the optimal conditions for your tubes and microwave.

1. Place 40µl DI water in the same type of tube that you will be using for GenePurgeDirect® treatment.
2. Overlay each tube with mineral oil to prevent evaporation.
3. Close the tubes, place in microwave safe rack (polyethylene or propylene) and heat on high for 5 minutes.
4. If any caps pop or tubes distort in any manner, then place a separate beaker in the microwave with 150ml of room temperature DI water and repeat the above 3 steps, the beaker of water serves as a heat ballast.
5. If tubes open or distort, reduce the power by 10% increments and increase time by 1-minute increments repeating step 4 until tubes no longer open or distort.

Note: Make sure the racks used in this procedure are MICROWAVE SAFE!

II. Microwave Protocol

1. Perform microwave procedure above for time and power conditions
1. Place 1µl of specimen with 20µl of GenePurgeDirect® into either a 0.5ml PCR tube or 1.5ml tube.
2. Vortex the tubes containing specimen and GenePurgeDirect® for ~10 seconds.
3. Overlay with mineral oil to prevent samples from evaporating.
4. Place the closed tubes in a microwave safe polyethylene or propylene rack. Make sure that the lids are loosely closed. If lids are closed too tightly tubes could rupture.
5. Place the rack in a microwave oven and heat at maximum power setting (setting should be based on the microwave evaluation results) for 5-7 minutes. Typically, 5 minutes if wattage is 900 or higher and 7 minutes if wattage is 500.
6. Remove rack from microwave and centrifuge the tubes at 5000xg for 5 minutes. After centrifuging samples, remove supernatant and use as DNA template.
7. Perform the amplification reaction

General Tips:

1. For best performance make sure GenePurgeDirect® is mixed thoroughly; it may need to be shaken between samples.
2. The specimen and GenePurgeDirect® volumes may be adjusted, however use no less than 5ul of GenePurgeDirect®.
3. GenePurgeDirect® can be UV treated for 30 minutes prior to use if required.
4. The volume of GenePurgeDirect® used to accomplish lysis should be compensated for by deducting an equivalent volume of water from the components of the amplification reagents to maintain their appropriate concentrations in the final reaction volume.
5. If upon use of this product, the expected bands are not observed, then a magnesium titration (1.5mM – 4.0mM) should be performed. If this fails to produce the desired bands, then a reduction in annealing temperature by 5°C to 10°C should be employed in conjunction with the magnesium titration.
6. If the GenePurgeDirect® treated specimens need to be stored prior to performing the amplification reaction, store at 4°C or -20°C until the treated specimens can be amplified. Prior to amplification, stored specimens should be heated to 80°C and the amplification begun using an initial denaturing step of at least 2 minutes.

Store at 4°C. DO NOT FREEZE.

No license to use PCR is given or to be implied by this product or these instructions. A license to use the PCR Process for certain research and development activities accompanies the purchase of licensed suppliers' PCR reagents when used in conjunction with an authorized thermal cycler.

References:

Arons E, Sorbara L, Raffeld M, Stetler-Stevenson M, Steinberg SM, Liewehr DJ, Pastan I, Kreitman RJ. Characterization of T-cell repertoire in hairy cell leukemia patients before and after recombinant immunotoxin BL22 therapy. *Cancer Immunology, Immunotherapy*. 2006 Sep;55(9):1100-10.

Dewhirst FE, Fox JG, Mendes EN, Paster BJ, Gates CE, Kirkbride CA, Eaton KA. 'Flexispira rappini' strains represent at least 10 *Helicobacter* taxa. *International Journal of Systematic and Evolutionary Microbiology* 2000; 50: 1781–1787.

Gwaltney S, Oberst R. Comparison of an improved polymerase chain reaction protocol and the indirect hemagglutination assay in the detection of Eperythrozoon suis infection. *Journal of Veterinary Diagnostics*. 1994 July;6:321-325.

Nerurkar VR, Nguyen HT, Dashwood WM, Hoffmann PR, Yin C, Morens DM, Kaplan AH, Detels R, Yanagihara R. HIV type 1 subtype E in commercial sex workers and injection drug users in southern Vietnam. *AIDS Research Human Retroviruses*. 1996 Jun 10;12(9):841-3.

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India
Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com