

# TFE3 Break Apart FISH Probe Kit


## Introduction

The TFE3 Break Apart FISH Probe Kit is designed to detect rearrangements in the human *TFE3* gene located on chromosome band Xp11.23. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other *TFE3* aberrations such as deletions or amplifications. Rearrangements and abnormal expression of the *TFE3* gene – also known as *TFEA*, *RCCP2*, *RCCX1* or *bHLHe33* – have been observed in renal cell carcinoma subtypes, alveolar soft part carcinoma, and other tumor types.

Intended Use
To detect rearrangements in the human <i>TFE3</i> gene located on chromosome band Xp11.23.

Cont.	Color
LSP TFE3 5' FISH Probe LSP TFE3 3' FISH Probe	CytoOrange CytoGreen


## Probe Design



Cat. No.	Volume
CT-PAC013-10-OG	10 Tests (100 µL)

Signal Pattern Interpretation	
<u>Normal Pattern</u> 2OG*	<u>Abnormal Pattern</u> Other Patterns
*Overlapping orange and green signals can appear as yellow.	

- Henthorn PS, et al. *Genomics*. 11(2):374-8 (1991).
- Numakura K, et al. *Int J Clin Oncol*. 16(5):577-80 (2011).
- Malouf GG, et al. *Ann Oncol*. 21(9):1834-8 (2010).
- Aulmann S, et al. *Histopathology*. 50(7):881-6 (2007).
- Rao Q, et al. *Am J Surg Pathol*. 37(6):804-15 (2013).



CytoTest Inc.  
9430 Key West Ave., Suite 210  
Rockville, MD 20850, USA

\* CE IVD only available in certain countries. All other countries are either ASR or RUO. Please contact your local dealer or our headquarters for more information.

### India Contact:

### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)