

## PathoScreen<sup>®</sup> Kit for TSV

---

**Test Label:** Alkaline Phosphatase

**Test Format:** Triple Antibody Sandwich (TAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Monoclonal Mixture

### Intended Use:

---

This ELISA is a qualitative serological assay intended for detection of *Tobacco streak virus* (TSV) in plant foliage.

The test demonstrates cross-reactivity with *Strawberry necrotic shock virus* (SNSV), therefore, can also be used as a detection tool for SNSV.

### Includes:

- Antibody-coated 96-well microtiter plates
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- ECI buffer, 5X concentrate
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- General extract buffer, GEB
- [User Guide](#)



PathoScreen<sup>®</sup> Kit, Tobacco streak virus (TSV),PSA 25500/0096

PathoScreen<sup>®</sup> Kit, Tobacco streak virus (TSV),PSA 25500/0288

PathoScreen<sup>®</sup> Kit, Tobacco streak virus (TSV),PSA 25500/0480

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)