

Veratox for Gliadin



Veratox for Gliadin/Gluten is used for the quantitative analysis of prolamins (wheat gliadin, rye secalin, and barley hordein) in ingredients, clean-in-place rinses and finished food products intended to be prolamins (gluten) free.


Features and Benefits

- Accurate and easy-to-use
- Cost effective microwell format
- High specificity

Product Specifications

- Range of quantitation: 5–50 ppm
- Controls provided:
 - 0, 5, 10, 20 and 50 ppm gliadin
 - 0, 10, 20, 40 and 100 ppm gluten
- Testing time: 30 minutes
- Tests per kit: Up to 38

Results



Order Details

Product #	Product Description
8480	Veratox for Gliadin R5
8429	Allergen Extraction Kit – 20 samples
8482	Gliadin Extraction additive for dark chocolate, cocoa and tannin
8483	Gliadin Cocktail Solution for heat processed samples – 100 mL
8483B	Gliadin Cocktail Solution for heat processed samples – 500 mL

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com