

ELISA Reagent Set for PVY common strains

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal Detection Reagent: Monoclonal

Intended Use:

This test is used for detection of common strains of PVY in foliage and sprouts of potato. It can also detect PVY (c strain) in several other solanaceous crops.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- User Guide



Reagent Set for Potato virus Y -c strain (PVY -c), 96 testwells, SRA 20700/0096

Reagent Set for Potato virus Y -c strain (PVY -c), 500 testwells, SRA 20700/0500

Reagent Set for Potato virus Y -c strain (PVY -c), 1000 testwells, SRA 20700/1000

Reagent Set for Potato virus Y -c strain (PVY -c), 5000 testwells, SRA 20700/5000

India Contact:

Life Technologies (India) Pvt. Ltd.