

## ELISA Reagent Set for Soybean vein necrosis virus (SVNV)

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

### Intended Use:

Agdia's SVNV ELISA is used for detection of *Soybean vein necrosis virus (SVNV)* in soybean.

### Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for Soybean vein necrosis virus (SVNV), 96 testwells, SRA 43200/0096

Reagent Set for Soybean vein necrosis virus (SVNV), 500 testwells, SRA 43200/0500

Reagent Set for Soybean vein necrosis virus (SVNV), 1000 testwells, SRA 43200/1000

Reagent Set for Soybean vein necrosis virus (SVNV), 5000 testwells, SRA 43200/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)