

**List of Primary Monoclonal Antibodies; Host: Mouse; Reactivity: Human;**

**Application: Flow Cytometry (FC)**

Catalog No.	Description	Conjugation/Label	Clone	Isotype	Size	Concentration
100201	Anti-CD2 Antibody	UC	OKT11	IgG1,k	100µg	0.5mg/mL
100203	Anti-CD2 Antibody	UC	OKT11	IgG1,k	500µg	0.5mg/mL
100214	Anti-CD2 Antibody	FITC	OKT11	IgG1,k	25tests	100µg/mL
100215	Anti-CD2 Antibody	FITC	OKT11	IgG1,k	100tests	100µg/mL
100216	Anti-CD2 Antibody	FITC	OKT11	IgG1,k	200tests	100µg/mL
100224	Anti-CD2 Antibody	PE	OKT11	IgG1,k	25tests	50µg/mL
100225	Anti-CD2 Antibody	PE	OKT11	IgG1,k	100tests	50µg/mL
100226	Anti-CD2 Antibody	PE	OKT11	IgG1,k	200tests	50µg/mL
100234	Anti-CD2 Antibody	PerCP	OKT11	IgG1,k	25tests	50µg/mL
100235	Anti-CD2 Antibody	PerCP	OKT11	IgG1,k	100tests	50µg/mL
100236	Anti-CD2 Antibody	PerCP	OKT11	IgG1,k	200tests	50µg/mL
100244	Anti-CD2 Antibody	APC	OKT11	IgG1,k	25tests	50µg/mL
100245	Anti-CD2 Antibody	APC	OKT11	IgG1,k	100tests	50µg/mL
100246	Anti-CD2 Antibody	APC	OKT11	IgG1,k	200tests	50µg/mL
100251	Anti-CD2 Antibody	Biotin	OKT11	IgG1,k	100µg	0.5mg/mL
100301	Anti-CD3 Antibody	UC	OKT3	IgG2a, k	100µg	0.5mg/mL
100303	Anti-CD3 Antibody	UC	OKT3	IgG2a,k	500µg	0.5mg/mL
100314	Anti-CD3 Antibody	FITC	OKT3	IgG2a,k	25tests	200µg/mL
100315	Anti-CD3 Antibody	FITC	OKT3	IgG2a,k	100tests	200µg/mL
100316	Anti-CD3 Antibody	FITC	OKT3	IgG2a, k	200tests	200µg/mL
100324	Anti-CD3 Antibody	PE	OKT3	IgG2a,k	25tests	12µg/mL
100325	Anti-CD3 Antibody	PE	OKT3	IgG2a,k	100tests	12µg/mL
100326	Anti-CD3 Antibody	PE	OKT3	IgG2a,k	200tests	12µg/mL
100334	Anti-CD3 Antibody	PerCP	OKT3	IgG2a,k	25tests	50µg/mL
100335	Anti-CD3 Antibody	PerCP	OKT3	IgG2a,k	100tests	50µg/mL
100336	Anti-CD3 Antibody	PerCP	OKT3	IgG2a,k	200tests	50µg/mL
100344	Anti-CD3 Antibody	APC	OKT3	IgG2a,k	25tests	12µg/mL
100345	Anti-CD3 Antibody	APC	OKT3	IgG2a,k	100tests	12µg/mL
100346	Anti-CD3 Antibody	APC	OKT3	IgG2a,k	200tests	12µg/mL
100351	Anti-CD3 Antibody	Biotin	OKT3	IgG2a,k	100µg	0.5mg/mL
100401	Anti-CD4 Antibody	UC	OKT4	IgG2b, k	100µg	0.5mg/mL
100403	Anti-CD4 Antibody	UC	OKT4	IgG2b, k	500µg	0.5mg/mL
100414	Anti-CD4 Antibody	FITC	OKT4	IgG2b, k	25tests	50µg/mL
100415	Anti-CD4 Antibody	FITC	OKT4	IgG2b, k	100tests	50µg/mL
100416	Anti-CD4 Antibody	FITC	OKT4	IgG2b, k	200tests	50µg/mL
100424	Anti-CD4 Antibody	PE	OKT4	IgG2b, k	25tests	50µg/mL
100425	Anti-CD4 Antibody	PE	OKT4	IgG2b, k	100tests	50µg/mL
100426	Anti-CD4 Antibody	PE	OKT4	IgG2b, k	200tests	50µg/mL
100434	Anti-CD4 Antibody	PerCP	OKT4	IgG2b, k	25tests	50µg/mL
100435	Anti-CD4 Antibody	PerCP	OKT4	IgG2b, k	100tests	50µg/mL
100436	Anti-CD4 Antibody	PerCP	OKT4	IgG2b, k	200tests	50µg/mL
100444	Anti-CD4 Antibody	APC	OKT4	IgG2b, k	25tests	25µg/mL
100445	Anti-CD4 Antibody	APC	OKT4	IgG2b, k	100tests	25µg/mL
100446	Anti-CD4 Antibody	APC	OKT4	IgG2b, k	200tests	25µg/mL
100451	Anti-CD4 Antibody	Biotin	OKT4	IgG2b, k	100µg	0.5mg/mL
100501	Anti-CD5 Antibody	UC	OKT1	IgG1,k	100µg	0.5mg/mL
100503	Anti-CD5 Antibody	UC	OKT1	IgG1,k	500µg	0.5mg/mL
100514	Anti-CD5 Antibody	FITC	OKT1	IgG1,k	25tests	100µg/mL
100515	Anti-CD5 Antibody	FITC	OKT1	IgG1,k	100tests	100µg/mL
100516	Anti-CD5 Antibody	FITC	OKT1	IgG1,k	200tests	100µg/mL
100524	Anti-CD5 Antibody	PE	OKT1	IgG1,k	25tests	50µg/mL
100525	Anti-CD5 Antibody	PE	OKT1	IgG1,k	100tests	50µg/mL
100526	Anti-CD5 Antibody	PE	OKT1	IgG1,k	200tests	50µg/mL
100534	Anti-CD5 Antibody	PerCP	OKT1	IgG1,k	25tests	50µg/mL
100535	Anti-CD5 Antibody	PerCP	OKT1	IgG1,k	100tests	50µg/mL

**List of Primary Monoclonal Antibodies; Host: Mouse; Reactivity: Human;**

**Application: Flow Cytometry (FC)**

100536	Anti-CD5 Antibody	PerCP	OKT1	IgG1,k	200tests	50µg/mL
100544	Anti-CD5 Antibody	APC	OKT1	IgG1,k	25tests	50µg/mL
100545	Anti-CD5 Antibody	APC	OKT1	IgG1,k	100tests	50µg/mL
100546	Anti-CD5 Antibody	APC	OKT1	IgG1,k	200tests	50µg/mL
100551	Anti-CD5 Antibody	Biotin	OKT1	IgG1,k	100µg	0.5mg/mL
100601	Anti-CD8 Antibody	UC	OKT8	IgG2a,k	100µg	0.5mg/mL
100603	Anti-CD8 Antibody	UC	OKT8	IgG2a,k	500µg	0.5mg/mL
100614	Anti-CD8 Antibody	FITC	OKT8	IgG2a,k	25tests	100µg/mL
100615	Anti-CD8 Antibody	FITC	OKT8	IgG2a,k	100tests	100µg/mL
100616	Anti-CD8 Antibody	FITC	OKT8	IgG2a,k	200tests	100µg/mL
100624	Anti-CD8 Antibody	PE	OKT8	IgG2a,k	25tests	50µg/mL
100625	Anti-CD8 Antibody	PE	OKT8	IgG2a,k	100tests	50µg/mL
100626	Anti-CD8 Antibody	PE	OKT8	IgG2a,k	200tests	50µg/mL
100634	Anti-CD8 Antibody	PerCP	OKT8	IgG2a,k	25tests	25µg/mL
100635	Anti-CD8 Antibody	PerCP	OKT8	IgG2a,k	100tests	25µg/mL
100636	Anti-CD8 Antibody	PerCP	OKT8	IgG2a,k	200tests	25µg/mL
100644	Anti-CD8 Antibody	APC	OKT8	IgG2a,k	25tests	12µg/mL
100645	Anti-CD8 Antibody	APC	OKT8	IgG2a,k	100tests	12µg/mL
100646	Anti-CD8 Antibody	APC	OKT8	IgG2a,k	200tests	12µg/mL
100651	Anti-CD8 Antibody	Biotin	OKT8	IgG2a,k	100µg	0.5mg/mL
100701	Anti-CD11b Antibody	UC	OKM1	IgG2b,k	100µg	0.5mg/mL
100703	Anti-CD11b Antibody	UC	OKM1	IgG2b,k	500µg	0.5mg/mL
100714	Anti-CD11b Antibody	FITC	OKM1	IgG2b,k	25tests	50µg/mL
100715	Anti-CD11b Antibody	FITC	OKM1	IgG2b,k	100tests	50µg/mL
100716	Anti-CD11b Antibody	FITC	OKM1	IgG2b,k	200tests	50µg/mL
100724	Anti-CD11b Antibody	PE	OKM1	IgG2b,k	25tests	50µg/mL
100725	Anti-CD11b Antibody	PE	OKM1	IgG2b,k	100tests	50µg/mL
100726	Anti-CD11b Antibody	PE	OKM1	IgG2b,k	200tests	50µg/mL
100734	Anti-CD11b Antibody	PerCP	OKM1	IgG2b,k	25tests	50µg/mL
100735	Anti-CD11b Antibody	PerCP	OKM1	IgG2b,k	100tests	50µg/mL
100736	Anti-CD11b Antibody	PerCP	OKM1	IgG2b,k	200tests	50µg/mL
100744	Anti-CD11b Antibody	APC	OKM1	IgG2b,k	25tests	35µg/mL
100745	Anti-CD11b Antibody	APC	OKM1	IgG2b,k	100tests	35µg/mL
100746	Anti-CD11b Antibody	APC	OKM1	IgG2b,k	200tests	35µg/mL
100751	Anti-CD11b Antibody	Biotin	OKM1	IgG2b,k	100µg	0.5mg/mL
100801	Anti-CD38 Antibody	UC	OKT10	IgG1,k	100µg	0.5mg/mL
100803	Anti-CD38 Antibody	UC	OKT10	IgG1,k	500µg	0.5mg/mL
100814	Anti-CD38 Antibody	FITC	OKT10	IgG1,k	25tests	200µg/mL
100815	Anti-CD38 Antibody	FITC	OKT10	IgG1,k	100tests	200µg/mL
100816	Anti-CD38 Antibody	FITC	OKT10	IgG1,k	200tests	200µg/mL
100824	Anti-CD38 Antibody	PE	OKT10	IgG1,k	25tests	50µg/mL
100825	Anti-CD38 Antibody	PE	OKT10	IgG1,k	100tests	50µg/mL
100826	Anti-CD38 Antibody	PE	OKT10	IgG1,k	200tests	50µg/mL
100851	Anti-CD38 Antibody	Biotin	OKT10	IgG1,k	100µg	0.5mg/mL
101601	Anti-CD45 Antibody	UC	F10-89-4	IgG2a,k	100µg	0.5mg/mL
101603	Anti-CD45 Antibody	UC	F10-89-4	IgG2a,k	500µg	0.5mg/mL
101614	Anti-CD45 Antibody	FITC	F10-89-4	IgG2a,k	25tests	50µg/mL
101615	Anti-CD45 Antibody	FITC	F10-89-4	IgG2a,k	100tests	50µg/mL
101616	Anti-CD45 Antibody	FITC	F10-89-4	IgG2a,k	200tests	50µg/mL
101624	Anti-CD45 Antibody	PE	F10-89-4	IgG2a,k	25tests	50µg/mL
101625	Anti-CD45 Antibody	PE	F10-89-4	IgG2a,k	100tests	50µg/mL
101626	Anti-CD45 Antibody	PE	F10-89-4	IgG2a,k	200tests	50µg/mL
101634	Anti-CD45 Antibody	PerCP	F10-89-4	IgG2a,k	25tests	50µg/mL
101635	Anti-CD45 Antibody	PerCP	F10-89-4	IgG2a,k	100tests	50µg/mL
101636	Anti-CD45 Antibody	PerCP	F10-89-4	IgG2a,k	200tests	50µg/mL
101644	Anti-CD45 Antibody	APC	F10-89-4	IgG2a,k	25tests	12µg/mL
101645	Anti-CD45 Antibody	APC	F10-89-4	IgG2a,k	100tests	12µg/mL
101646	Anti-CD45 Antibody	APC	F10-89-4	IgG2a,k	200tests	12µg/mL
101651	Anti-CD45 Antibody	Biotin	F10-89-4	IgG2a,k	100µg	0.5mg/mL

**List of Primary Monoclonal Antibodies; Host: Mouse; Reactivity: Human;**

**Application: Flow Cytometry (FC)**

101801	Anti-CD90 Antibody	UC	F15-42-1	IgG1,k	100µg	0.5mg/mL
101803	Anti-CD90 Antibody	UC	F15-42-1	IgG1,k	500µg	0.5mg/mL
101814	Anti-CD90 Antibody	FITC	F15-42-1	IgG1,k	25tests	200µg/mL
101815	Anti-CD90 Antibody	FITC	F15-42-1	IgG1,k	100tests	200µg/mL
101816	Anti-CD90 Antibody	FITC	F15-42-1	IgG1,k	200tests	200µg/mL
101824	Anti-CD90 Antibody	PE	F15-42-1	IgG1,k	25tests	50µg/mL
101825	Anti-CD90 Antibody	PE	F15-42-1	IgG1,k	100tests	50µg/mL
101826	Anti-CD90 Antibody	PE	F15-42-1	IgG1,k	200tests	50µg/mL
101844	Anti-CD90 Antibody	APC	F15-42-1	IgG1,k	25tests	50µg/mL
101845	Anti-CD90 Antibody	APC	F15-42-1	IgG1,k	100tests	50µg/mL
101846	Anti-CD90 Antibody	APC	F15-42-1	IgG1,k	200tests	50µg/mL
101851	Anti-CD90 Antibody	Biotin	F15-42-1	IgG1,k	100µg	0.5mg/mL
102301	Isotype Control IgG1,k	UC	MOPC-31C	IgG1,k	100µg	0.5mg/mL
102303	Isotype Control IgG1,k	UC	MOPC-31C	IgG1,k	500µg	0.5mg/mL
102314	Isotype Control IgG1,k	FITC	MOPC-31C	IgG1,k	25tests	100µg/mL
102315	Isotype Control IgG1,k	FITC	MOPC-31C	IgG1,k	100tests	100µg/mL
102316	Isotype Control IgG1,k	FITC	MOPC-31C	IgG1,k	200tests	100µg/mL
102324	Isotype Control IgG1,k	PE	MOPC-31C	IgG1,k	25tests	50µg/mL
102325	Isotype Control IgG1,k	PE	MOPC-31C	IgG1,k	100tests	50µg/mL
102326	Isotype Control IgG1,k	PE	MOPC-31C	IgG1,k	200tests	50µg/mL
102335	Isotype Control IgG1,k	PerCP	MOPC-31C	IgG1,k	100tests	50µg/mL
102336	Isotype Control IgG1,k	PerCP	MOPC-31C	IgG1,k	200tests	50µg/mL
102344	Isotype Control IgG1,k	APC	MOPC-31C	IgG1,k	25tests	50µg/mL
102345	Isotype Control IgG1,k	APC	MOPC-31C	IgG1,k	100tests	50µg/mL
102346	Isotype Control IgG1,k	APC	MOPC-31C	IgG1,k	200tests	50µg/mL
102351	Isotype Control IgG1,k	Biotin	MOPC-31C	IgG1,k	100µg	0.5mg/mL
102401	Isotype Control IgG2b,k	UC	MPC-11	IgG2b,k	100µg	0.5mg/mL
102403	Isotype Control IgG2b,k	UC	MPC-11	IgG2b,k	500µg	0.5mg/mL
102414	Isotype Control IgG2b,k	FITC	MPC-11	IgG2b,k	25tests	100µg/mL
102415	Isotype Control IgG2b,k	FITC	MPC-11	IgG2b,k	100tests	100µg/mL
102416	Isotype Control IgG2b,k	FITC	MPC-11	IgG2b,k	200tests	100µg/mL
102424	Isotype Control IgG2b,k	PE	MPC-11	IgG2b,k	25tests	50µg/mL
102425	Isotype Control IgG2b,k	PE	MPC-11	IgG2b,k	100tests	50µg/mL
102426	Isotype Control IgG2b,k	PE	MPC-11	IgG2b,k	200tests	50µg/mL
102435	Isotype Control IgG2b,k	PerCP	MPC-11	IgG2b,k	100tests	50µg/mL
102436	Isotype Control IgG2b,k	PerCP	MPC-11	IgG2b,k	200tests	50µg/mL
102444	Isotype Control IgG2b,k	APC	MPC-11	IgG2b,k	25tests	50µg/mL
102445	Isotype Control IgG2b,k	APC	MPC-11	IgG2b,k	100tests	50µg/mL
102446	Isotype Control IgG2b,k	APC	MPC-11	IgG2b,k	200tests	50µg/mL
102451	Isotype Control IgG2b,k	Biotin	MPC-11	IgG2b,k	100µg	0.5mg/mL
102601	Anti-CD25 Antibody	UC	7G7B6	IgG2a,k	100µg	0.5mg/mL
102603	Anti-CD25 Antibody	UC	7G7B6	IgG2a,k	500µg	0.5mg/mL
102614	Anti-CD25 Antibody	FITC	7G7B6	IgG2a,k	25tests	100µg/mL
102615	Anti-CD25 Antibody	FITC	7G7B6	IgG2a,k	100tests	100µg/mL
102616	Anti-CD25 Antibody	FITC	7G7B6	IgG2a,k	200tests	100µg/mL
102624	Anti-CD25 Antibody	PE	7G7B6	IgG2a,k	25tests	50µg/mL
102625	Anti-CD25 Antibody	PE	7G7B6	IgG2a,k	100tests	50µg/mL
102626	Anti-CD25 Antibody	PE	7G7B6	IgG2a,k	200tests	50µg/mL
102651	Anti-CD25 Antibody	Biotin	7G7B6	IgG2a,k	100µg	0.5mg/mL
102701	Anti-CD32 Antibody	UC	IV.3	IgG2b,k	100µg	0.5mg/mL
102703	Anti-CD32 Antibody	UC	IV.3	IgG2b,k	500µg	0.5mg/mL
102714	Anti-CD32 Antibody	FITC	IV.3	IgG2b,k	25tests	25µg/mL
102715	Anti-CD32 Antibody	FITC	IV.3	IgG2b,k	100tests	25µg/mL
102716	Anti-CD32 Antibody	FITC	IV.3	IgG2b,k	200tests	25µg/mL
102724	Anti-CD32 Antibody	PE	IV.3	IgG2b,k	25tests	50µg/mL
102725	Anti-CD32 Antibody	PE	IV.3	IgG2b,k	100tests	50µg/mL
102726	Anti-CD32 Antibody	PE	IV.3	IgG2b,k	200tests	50µg/mL
102744	Anti-CD32 Antibody	APC	IV.3	IgG2b,k	25tests	25µg/mL
102745	Anti-CD32 Antibody	APC	IV.3	IgG2b,k	100tests	25µg/mL

**List of Primary Monoclonal Antibodies; Host: Mouse; Reactivity: Human;**

**Application: Flow Cytometry (FC)**

102746	Anti-CD32 Antibody	APC	IV.3	IgG2b,k	200tests	25µg/mL
102751	Anti-CD32 Antibody	Biotin	IV.3	IgG2b,k	100µg	0.5mg/mL
102801	Anti-CD62L Antibody	UC	DREG56	IgG1,k	100µg	0.5mg/mL
102803	Anti-CD62L Antibody	UC	DREG56	IgG1,k	500µg	0.5mg/mL
102814	Anti-CD62L Antibody	FITC	DREG56	IgG1,k	25tests	25µg/mL
102815	Anti-CD62L Antibody	FITC	DREG56	IgG1,k	100tests	25µg/mL
102816	Anti-CD62L Antibody	FITC	DREG56	IgG1,k	200tests	25µg/mL
102824	Anti-CD62L Antibody	PE	DREG56	IgG1,k	25tests	25µg/mL
102825	Anti-CD62L Antibody	PE	DREG56	IgG1,k	100tests	25µg/mL
102826	Anti-CD62L Antibody	PE	DREG56	IgG1,k	200tests	25µg/mL
102844	Anti-CD62L Antibody	APC	DREG56	IgG1,k	25tests	12µg/mL
102845	Anti-CD62L Antibody	APC	DREG56	IgG1,k	100tests	12µg/mL
102846	Anti-CD62L Antibody	APC	DREG56	IgG1,k	200tests	12µg/mL
102851	Anti-CD62L Antibody	Biotin	DREG56	IgG1,k	100µg	0.5mg/mL
102901	Anti-CD19 Antibody	UC	4G7	IgG1, k	100µg	0.5mg/mL
102903	Anti-CD19 Antibody	UC	4G7	IgG1, k	500µg	0.5mg/mL
102914	Anti-CD19 Antibody	FITC	4G7	IgG1, k	25tests	100µg/mL
102915	Anti-CD19 Antibody	FITC	4G7	IgG1, k	100tests	100µg/mL
102916	Anti-CD19 Antibody	FITC	4G7	IgG1, k	200tests	100µg/mL
102924	Anti-CD19 Antibody	PE	4G7	IgG1, k	25tests	25µg/mL
102925	Anti-CD19 Antibody	PE	4G7	IgG1, k	100tests	25µg/mL
102926	Anti-CD19 Antibody	PE	4G7	IgG1, k	200tests	25µg/mL
102934	Anti-CD19 Antibody	PerCP	4G7	IgG1, k	25tests	50µg/mL
102935	Anti-CD19 Antibody	PerCP	4G7	IgG1, k	100tests	50µg/mL
102936	Anti-CD19 Antibody	PerCP	4G7	IgG1, k	200tests	50µg/mL
102944	Anti-CD19 Antibody	APC	4G7	IgG1, k	25tests	25µg/mL
102945	Anti-CD19 Antibody	APC	4G7	IgG1, k	100tests	25µg/mL
102946	Anti-CD19 Antibody	APC	4G7	IgG1, k	200tests	25µg/mL
102951	Anti-CD19 Antibody	Biotin	4G7	IgG1, k	100µg	0.5mg/mL
103001	Anti-CD7 Antibody	UC	4H9	IgG2a,k	100µg	0.5mg/mL
103003	Anti-CD7 Antibody	UC	4H9	IgG2a,k	500µg	0.5mg/mL
103014	Anti-CD7 Antibody	FITC	4H9	IgG2a,k	25tests	200µg/mL
103015	Anti-CD7 Antibody	FITC	4H9	IgG2a,k	100tests	200µg/mL
103016	Anti-CD7 Antibody	FITC	4H9	IgG2a,k	200tests	200µg/mL
103024	Anti-CD7 Antibody	PE	4H9	IgG2a,k	25tests	50µg/mL
103025	Anti-CD7 Antibody	PE	4H9	IgG2a,k	100tests	50µg/mL
103026	Anti-CD7 Antibody	PE	4H9	IgG2a,k	200tests	50µg/mL
103044	Anti-CD7 Antibody	APC	4H9	IgG2a,k	25tests	50µg/mL
103045	Anti-CD7 Antibody	APC	4H9	IgG2a,k	100tests	50µg/mL
103046	Anti-CD7 Antibody	APC	4H9	IgG2a,k	200tests	50µg/mL
103051	Anti-CD7 Antibody	Biotin	4H9	IgG2a,k	100µg	0.5mg/mL
103101	Anti-CD71 Antibody	UC	OKT9	IgG1,k	100µg	0.5mg/mL
103103	Anti-CD71 Antibody	UC	OKT9	IgG1,k	500µg	0.5mg/mL
103114	Anti-CD71 Antibody	FITC	OKT9	IgG1,k	25tests	100µg/mL
103115	Anti-CD71 Antibody	FITC	OKT9	IgG1,k	100tests	100µg/mL
103116	Anti-CD71 Antibody	FITC	OKT9	IgG1,k	200tests	100µg/mL
103124	Anti-CD71 Antibody	PE	OKT9	IgG1,k	25tests	25µg/mL
103125	Anti-CD71 Antibody	PE	OKT9	IgG1,k	100tests	25µg/mL
103126	Anti-CD71 Antibody	PE	OKT9	IgG1,k	200tests	25µg/mL
103144	Anti-CD71 Antibody	APC	OKT9	IgG1,k	25tests	12µg/mL
103145	Anti-CD71 Antibody	APC	OKT9	IgG1,k	100tests	12µg/mL
103146	Anti-CD71 Antibody	APC	OKT9	IgG1,k	200tests	12µg/mL
103151	Anti-CD71 Antibody	Biotin	OKT9	IgG1,k	100µg	0.5mg/mL
103201	Anti-HLA-DR Antibody	UC	L243	IgG2a,k	100µg	0.5mg/mL
103203	Anti-HLA-DR Antibody	UC	L243	IgG2a,k	500µg	0.5mg/mL
103214	Anti-HLA-DR Antibody	FITC	L243	IgG2a,k	25tests	100µg/mL
103215	Anti-HLA-DR Antibody	FITC	L243	IgG2a,k	100tests	100µg/mL
103216	Anti-HLA-DR Antibody	FITC	L243	IgG2a,k	200tests	100µg/mL
103224	Anti-HLA-DR Antibody	PE	L243	IgG2a,k	25tests	50µg/mL

**List of Primary Monoclonal Antibodies; Host: Mouse; Reactivity: Human;**

**Application: Flow Cytometry (FC)**

103225	Anti-HLA-DR Antibody	PE	L243	IgG2a,k	100tests	50µg/mL
103226	Anti-HLA-DR Antibody	PE	L243	IgG2a,k	200tests	50µg/mL
103234	Anti-HLA-DR Antibody	PerCP	L243	IgG2a,k	25tests	50µg/mL
103235	Anti-HLA-DR Antibody	PerCP	L243	IgG2a,k	100tests	50µg/mL
103236	Anti-HLA-DR Antibody	PerCP	L243	IgG2a,k	200tests	50µg/mL
103244	Anti-HLA-DR Antibody	APC	L243	IgG2a,k	25tests	12µg/mL
103245	Anti-HLA-DR Antibody	APC	L243	IgG2a,k	100tests	12µg/mL
103246	Anti-HLA-DR Antibody	APC	L243	IgG2a,k	200tests	12µg/mL
103251	Anti-HLA-DR Antibody	Biotin	L243	IgG2a,k	100µg	0.5mg/mL
103301	Anti-CD57 Antibody	UC	HNK-1	IgM	100µg	0.5mg/mL
103301	Anti-CD57 Antibody	UC	HNK-1	IgM	500µg	0.5mg/mL
103324	Anti-CD57 Antibody	PE	HNK-1	IgM	25tests	25µg/mL
103325	Anti-CD57 Antibody	PE	HNK-1	IgM	100tests	25µg/mL
103326	Anti-CD57 Antibody	PE	HNK-1	IgM	200tests	25µg/mL
103351	Anti-CD57 Antibody	Biotin	HNK-1	IgM	100µg	0.5mg/mL
103401	Anti-CD14 Antibody	UC	26ic	IgG2b, k	100µg	0.5mg/mL
103403	Anti-CD14 Antibody	UC	26ic	IgG2b, k	500µg	0.5mg/mL
103414	Anti-CD14 Antibody	FITC	26ic	IgG2b, k	25tests	25µg/mL
103415	Anti-CD14 Antibody	FITC	26ic	IgG2b, k	100tests	25µg/mL
103416	Anti-CD14 Antibody	FITC	26ic	IgG2b, k	200tests	25µg/mL
103424	Anti-CD14 Antibody	PE	26ic	IgG2b, k	25tests	25µg/mL
103425	Anti-CD14 Antibody	PE	26ic	IgG2b, k	100tests	25µg/mL
103426	Anti-CD14 Antibody	PE	26ic	IgG2b, k	200tests	25µg/mL
103434	Anti-CD14 Antibody	PerCP	26ic	IgG2b, k	25tests	12µg/mL
103435	Anti-CD14 Antibody	PerCP	26ic	IgG2b, k	100tests	12µg/mL
103436	Anti-CD14 Antibody	PerCP	26ic	IgG2b, k	200tests	12µg/mL
103444	Anti-CD14 Antibody	APC	26ic	IgG2b, k	25tests	12µg/mL
103445	Anti-CD14 Antibody	APC	26ic	IgG2b, k	100tests	12µg/mL
103446	Anti-CD14 Antibody	APC	26ic	IgG2b, k	200tests	12µg/mL
103451	Anti-CD14 Antibody	Biotin	26ic	IgG2b, k	100µg	0.5mg/mL
103524	Anti-CD117 Antibody	PE	BA7.3C.9	IgG2a,k	25tests	25µg/mL
103525	Anti-CD117 Antibody	PE	BA7.3C.9	IgG2a,k	100tests	25µg/mL
103526	Anti-CD117 Antibody	PE	BA7.3C.9	IgG2a,k	200tests	25µg/mL
103601	Isotype Control IgG2a,k	UC	RPC5.4	IgG2a,k	100µg	0.5mg/mL
103603	Isotype Control IgG2a,k	UC	RPC5.4	IgG2a,k	500µg	0.5mg/mL
103614	Isotype Control IgG2a,k	FITC	RPC5.4	IgG2a,k	25tests	100µg/mL
103615	Isotype Control IgG2a,k	FITC	RPC5.4	IgG2a,k	100tests	100µg/mL
103616	Isotype Control IgG2a,k	FITC	RPC5.4	IgG2a,k	200tests	100µg/mL
103624	Isotype Control IgG2a,k	PE	RPC5.4	IgG2a,k	25tests	50µg/mL
103625	Isotype Control IgG2a,k	PE	RPC5.4	IgG2a,k	100tests	50µg/mL
103626	Isotype Control IgG2a,k	PE	RPC5.4	IgG2a,k	200tests	50µg/mL
103635	Isotype Control IgG2a,k	PerCP	RPC5.4	IgG2a,k	100tests	50µg/mL
103636	Isotype Control IgG2a,k	PerCP	RPC5.4	IgG2a,k	200tests	50µg/mL
103644	Isotype Control IgG2a,k	APC	RPC5.4	IgG2a,k	25tests	50µg/mL
103645	Isotype Control IgG2a,k	APC	RPC5.4	IgG2a,k	100tests	50µg/mL
103646	Isotype Control IgG2a,k	APC	RPC5.4	IgG2a,k	200tests	50µg/mL
103651	Isotype Control IgG2a,k	Biotin	RPC5.4	IgG2a,k	100µg	0.5mg/mL
103701	Anti- CD20 Antibody	UC	2H7	IgG2b, k	100µg	0.5mg/mL
103703	Anti- CD20 Antibody	UC	2H7	IgG2b, k	500µg	0.5mg/mL
103714	Anti- CD20 Antibody	FITC	2H7	IgG2b, k	25tests	100µg/mL
103715	Anti- CD20 Antibody	FITC	2H7	IgG2b, k	100tests	100µg/mL
103716	Anti-CD20 Antibody	FITC	2H7	IgG2b, k	200tests	100µg/mL
103724	Anti- CD20 Antibody	PE	2H7	IgG2b, k	25tests	50µg/mL
103725	Anti- CD20 Antibody	PE	2H7	IgG2b, k	100tests	50µg/mL
103726	Anti-CD20 Antibody	PE	2H7	IgG2b, k	200tests	50µg/mL
103734	Anti-CD20 Antibody	PerCP	2H7	IgG2b, k	25tests	50µg/mL
103735	Anti-CD20 Antibody	PerCP	2H7	IgG2b, k	100tests	50µg/mL
103736	Anti-CD20 Antibody	PerCP	2H7	IgG2b, k	200tests	50µg/mL
103744	Anti-CD20 Antibody	APC	2H7	IgG2b, k	25tests	50µg/mL

**List of Primary Monoclonal Antibodies; Host: Mouse; Reactivity: Human;**

**Application: Flow Cytometry (FC)**

103745	Anti-CD20 Antibody	APC	2H7	IgG2b, k	100tests	50µg/mL
103746	Anti-CD20 Antibody	APC	2H7	IgG2b, k	200tests	50µg/mL
103751	Anti-CD20 Antibody	Biotin	2H7	IgG2b, k	100µg	0.5mg/mL
102101	Anti-p53 Antibody for ICC/WB	UC	PAB122	IgG2b, k	100µg	0.5mg/mL
100101	Anti-EAPII Antibody for ICC/WB	UC	42c	IgG2b, k	100µg	0.5mg/mL
102201	Anti-p53 WT Antibody for ICC	UC	PAB1620	IgG2b, k	100µg	0.5mg/mL
104101	Anti-CD55 Antibody	UC	F4-29D9	IgG1,k	100µg	0.5mg/mL
104103	Anti-CD55 Antibody	UC	F4-29D9	IgG1,k	500µg	0.5mg/mL
104114	Anti-CD55 Antibody	FITC	F4-29D9	IgG1,k	25tests	
104115	Anti-CD55 Antibody	FITC	F4-29D9	IgG1,k	100tests	
104116	Anti-CD55 Antibody	FITC	F4-29D9	IgG1,k	200tests	
104124	Anti-CD55 Antibody	PE	F4-29D9	IgG1,k	25tests	
104125	Anti-CD55 Antibody	PE	F4-29D9	IgG1,k	100tests	
104126	Anti-CD55 Antibody	PE	F4-29D9	IgG1,k	200tests	
104001	Anti-CD33 Antibody	UC	MW53	IgG1,k	100µg	0.5mg/mL
104001	Anti-CD33 Antibody	UC	MW54	IgG1,k	500µg	0.5mg/mL
104024	Anti-CD33 Antibody	FITC	MW55	IgG1,k	25tests	
104025	Anti-CD33 Antibody	FITC	MW56	IgG1,k	100tests	
104026	Anti-CD33 Antibody	FITC	MW57	IgG1,k	200tests	
104301	Anti-CD235a Antibody	UC	GYPA280	IgG1,k	100µg	0.5mg/mL
104303	Anti-CD235a Antibody	UC	GYPA281	IgG1,k	500µg	0.5mg/mL
104314	Anti-CD235a Antibody	FITC	GYPA282	IgG1,k	25µg	0.5mg/mL
104315	Anti-CD235a Antibody	FITC	GYPA283	IgG1,k	100µg	0.5mg/mL
104324	Anti-CD235a Antibody	PE	GYPA284	IgG1,k	25µg	0.05mg/ml
104325	Anti-CD235a Antibody	PE	GYPA285	IgG1,k	100µg	0.05mg/ml
104201	Anti-CD59 Antibody	UC	DRA-10G	IgG1,k	100µg	0.5mg/mL
104203	Anti-CD59 Antibody	UC	DRA-10G	IgG1,k	500µg	0.5mg/mL

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, [www.atzlabs.com](http://www.atzlabs.com) ; [www.lifetechindia.com](http://www.lifetechindia.com)