


## List of ELISA Kits and Reagents for Food Safety and Drug Residues

ADI has highly specialized Drug residue testing ELISA kits for animal, human, or aquatic samples.

Product Name	Cat#	# tests Assay/ Range (ppb)	Samples	Detection Limit (ppb)	Sample Processing Time	Assay time
Clenbuterol ELISA Test Kit	DE-100010	96 tests 0.1-8.1 ppb	Pork/Liver	0.1	50 min	80 min
			Serum	0.1	10 min	80 min
			Urine	0.1	10 min	80 min
			Feed	50	30 min	80 min
			Meat	0.2	20 min	80 min
Ractopamine ELISA Test Kit	DE-100020	96 tests 0.1-8.1 ppb	Liver	0.4	20 min	80 min
			Urine	0.2	10 min	80 min
			Feed	15	30 min	50 min
Salbutamal ELISA Test Kit	DE-100030	96 tests 0.5-40.5 ppb	Urine	0.5	10 min	80 min
			Tissue	0.5	20 min	80 min
			Feed	50	30 min	80 min
Chloramphenicol ELISA Test Kit	DE-100040	96 tests 0.05-4.05 ppb	Animal tissue	0.05	20 min	80 min
			Aquatic	0.05	20 min	80 min
			Honey	0.15	20 min	80 min
			Intestine	0.1	20 min	80 min
			Urine	0.1	150 min	80 min
			Egg	0.1	30 min	80 min
			Milk	0.05	30 min	80 min
			Serum	0.1	15 min	80 min
Florfenicol ELISA Test Kit	DE-100050	96 tests 0.5-40.5 ppb	Animal tissue	0.5	30 min	80 min
			Aquatic	0.5	30 min	80 min
			Honey	1.5	30 min	80 min
			Fish, Shrimp	0.3	120 min	80 min
Nitrofurans (AMOZ) ELISA Test Kit	DE-100060	96 tests 0.1-8.1 ppb	Honey	0.2	120 min	80 min
			Chicken/Liver	0.2	120 min	80 min
Nitrofurans (AHD) ELISA Test Kit	DE-100070	96 tests 0.05-4.05 ppb	Fish, Shrimp	0.3	16 h	80 min
			Honey	0.2	16 h	80 min
			Chicken/Liver	0.2	16 h	80 min
Nitrofurans (AOZ) ELISA Test Kit	DE-100080	96 tests 0.05-4.05 ppb	Fish, Shrimp	0.2	16 h	80 min
			Honey	0.1	16 h	80 min
			Chicken/Liver	0.1	16 h	80 min
Sulfonamides residues (SAs) ELISA Test Kit	DE-100090	96 tests 1-81 ppb	Chicken/Liver	1	20 min	90 min
			Pork/Liver	1	20 min	90 min
			Honey/Egg	1	10 min	90 min
			Serum/Urine	4	10 min	90 min
			Milk	20	10 min	90 min
Sulfamethazine (SM2) ELISA Test Kit	DE-100100	96 tests 1-81 ppb	Chicken/Liver	1	20 min	90 min
			Pork/Liver	1	20 min	90 min
			Honey/Egg	1	10 min	90 min
			Serum/Urine	4	10 min	90 min
			Milk	20	10 min	90 min
Sulfamethoxydiazine (SMD) ELISA Test Kit	DE-100110	96 tests 1-81 ppb	Chicken/Liver	1	20 min	90 min
			Pork/Liver	1	20 min	90 min
			Honey/Egg	1	10 min	90 min
			Serum/Urine	4	10 min	90 min
Sulfaquin-oxaline ELISA Test Kit	DE-100240	96 tests 1-81 ppb	Pork/Liver	1	20 min	90 min
			Honey/Egg	1	10 min	90 min
			Serum/Urine	4	10 min	90 min
			Milk	20	10 min	90 min



## List of ELISA Kits and Reagents for Food Safety and Drug Residues

ADI has highly specialized Drug residue testing ELISA kits for animal, human, or aquatic samples.

Product Name	Cat#	# tests Assay/ Range (ppb)	Samples	Detection Limit (ppb)	Sample Processing Time	Assay time
Quinolones (QNS) ELISA Test Kit	DE-100120	96 tests 1-81 ppb	Pork/Liver	3	50 min	90 min
			Chicken/Liver	3	50 min	90 min
			Shrimp, Fish	3	20 min	90 min
			Serum	3	10 min	90 min
			Honey	6	10 min	80 min
Enrofloxacin ELISA Test Kit	DE-100130	96 tests 0.25-20.25 ppb	Pork/Liver	1	50 min	80 min
			Chicken/ liver	1	50 min	80 min
			Shrimp, Fish	1	20 min	80 min
			Serum	1	10 min	80 min
			Honey	2	10 min	80 min
Ampicillin ELISA Test Kit	DE-100140	96 tests 0.1-8.1 ppb	Pork/Liver	2	15 min	80 min
			Chicken (meat)	2	15 min	80 min
			Duck (meat)	2	15 min	80 min
			Shrimp, Fish	2	15 min	80 min
			honey, milk	2	10 min	80 min
Benzyl penicillin ELISA Test Kit	DE-100150	96 tests 0.1-8.1 ppb	Pork/liver	2	15 min	80 min
			Chicken(meat, liver)	2	15 min	80 min
			Duck(meat, liver)	2	15 min	80 min
			Shrimp, Fish	2	15 min	80 min
			Honey, Milk	2	10 min	80 min
Tylosin ELISA Test Kit	DE-100160	96 tests 1.5-121.5 ppb	Meat	3	20 min	60 min
			Liver	3	20 min	60 min
			Honey/Egg	1.5	10 min	60 min
Trenbolone ELISA Test Kit	DE-100170	96 tests 0.1-8.1 ppb	Animal tissue	5	30 min	80 min
			Aquatic	5	30 min	80 min
			Urine	0.5	30 min	80 min
Diazepam ELISA Test Kit	DE-100180-	96 tests 0.1-8.1 ppb	Tissue	1	30 min	80 min
			Urine	1	30 min	80 min
			Feed	20	30 min	80 min
Diethylstilbestrol (DES) ELISA Test Kit	DE-100190	96 tests 0.1-8.1 ppb	Fish, Shrimp	0.2	20 min	80 min
			Liver, Meat	2	20 min	80 min
			Feed	20	20 min	80 min
			Urine	0.6	20 min	80 min
Gentamicin ELISA Test Kit	DE-100200	96 tests 0.1-8.1 ppb	Chicken/Liver	4	40 min	80 min
Streptomycin ELISA Test Kit	DE-100210	96 tests 0.1-25.6 ppb	Chicken/Liver	20	40 min	60 min
			Honey	10	40 min	60 min
			Milk	20	5 min	60 min
Nitrofurantoin (SEM) ELISA Test Kit	DE-100070	96 tests 0.1-8.1 ppb	Honey	0.2	16 h	80 min
			Fish, Shrimp	0.3	16 h	80 min
			Chicken/Liver	0.2	16 h	80 min
			Fish, Shrimp	0.3	16 h	80 min
			Chicken/Liver	0.2	16 h	80 min
Olaquinox ELISA Test Kit	DE-100230	96 tests 1-81 ppb	Tissue	1	40 min	80 min
			Feed	100	45 min	80 min

Rev. 130118A



India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com) Website: [www.atzlabs.com](http://www.atzlabs.com)