

BioKits BLG Assay Kit



BioKits BLG Assay Kit is used for the detection and quantification of β -lactoglobulin (BLG) in food products and environmental swabs by enzyme immunoassay.

Product Specifications

- Lower limit of detection: 2 ppm BLG
- Range of quantitation: 2.5–40 ppm
- Testing time: 120 minutes
- Wells per kit: 48

Results



Order Details

Product #Product Description

902061Y BioKits BLG Assay Kit

901042J BioKits Allergen Swabbing Kit – Contains 20 swabs, tubes and buffer

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com