

# BioKits Peanut Assay Kit



BioKits Peanut Assay Kit is used for the detection and quantification of peanut in foods and environmental swabs.

## Product Specifications

- Lower limit of detection: 0.1 ppm peanut
- Range of quantitation: 1–20 ppm
- Testing time: 75 minutes
- Wells per kit: 96

## Results



## Order Details

### Product #Product Description

902048Q BioKits Peanut Assay Kit

901042J BioKits Allergen Swabbing Kit – Contains 20 swabs, tubes and buffer

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com), [www.atzlabs.com](http://www.atzlabs.com) ; [www.lifetechindia.com](http://www.lifetechindia.com)