


TMPRSS2-ERG Fusion/Translocation FISH Probe Kit Cat No CT-PAC176

Prostate Cancer Diagnostic Kit

The TMPRSS2-ERG fusion FISH probe set is designed to detect rearrangements between the human *TMPRSS2* and *ERG* genes, both located on the long arm of chromosome 21 (21q22.3 and 21q22.2, respectively).*



The probe is labeled with CytoOrange and CytoGreen. CytoOrange is a fluorophore with an excitation peak at 551nm and emission peak at 575nm, giving it an orange fluorescent signal with an appropriate filter; CytoGreen has an excitation peak at 495nm and emission peak at 518nm, giving it a green fluorescent signal with an appropriate filter.

* CytoTest shall not be responsible for patent infringement or any other violation that may occur with the use of this product.