

Veratox for Oxytetracycline



Veratox® for Oxytetracycline is a competitive enzyme immunoassay for the quantitative analysis of Oxytetracycline in muscle.

Veratox for Oxytetracycline enables international and government regulatory agencies, food manufacturers and processors, as well as quality assurance organizations, to detect Oxytetracycline in muscle and to satisfy customer concerns about food safety. Veratox for Oxytetracycline is an accurate, rapid, cost-effective and easy-to-use diagnostic tool compared to instrument based methods like HPLC and GC.

Product Specifications

- Range of quantitation: 3.6-324 ppb
- Controls provided: 0, 0.15, 0.75, 1.5, 4.5, and 13.5 ng/mL Oxytetracycline
- Testing time: 90 minutes
- Test per kit: Up to 36 when run in duplicate

Results



Order Details

Product #	Product Description
DR081	Veratox for Oxytetracycline

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com