

PathoScreen[®] Kit for INSV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Impatiens necrotic spot virus* (INSV) in plant foliage with or without external disease symptoms.

The test is proven to detect over 18 INSV isolates from all over the world.

Includes:

- Antibody-coated 96-well microtiter plates
- Alkaline phosphatase enzyme conjugate
- ECI buffer, 5X concentrate
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- General extract buffer, GEB
- [User Guide](#)



PathoScreen[®] Kit, Impatiens necrotic spot virus (INSV), PSA 20500/0096

PathoScreen[®] Kit, Impatiens necrotic spot virus (INSV), PSA 20500/0288

PathoScreen[®] Kit, Impatiens necrotic spot virus (INSV), PSA 20500/0480

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com