

## PathoScreen<sup>®</sup> Kit for ArMV

---

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

### Intended Use:

---

This ELISA is a qualitative serological assay for detection of *Arabidopsis mosaic virus*(ArMV) in plant foliage.

The test is proven to detect over 33 ArMV isolates from all over the world.

The test has been evaluated by Agdia Testing Service Department in over 162 different plant hosts including ornamental, vegetable, and fruit species. Agdia ArMV ELISA test is used in our Testing Services Ornamental, Berry, Petunia, Hosta, and Rose crop screens.

### Includes:

- Antibody-coated 96-well microtiter plates
- Alkaline phosphatase enzyme conjugate
- RUB6 enzyme conjugate diluent, 1X ready to use
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- General extract buffer 3, GEB3
- [User Guide](#)



PathoScreen<sup>®</sup> Kit, Arabis mosaic virus (ArMV),PSA 23203/0096

PathoScreen<sup>®</sup> Kit, Arabis mosaic virus (ArMV),PSA 23203/0288

PathoScreen<sup>®</sup> Kit, Arabis mosaic virus (ArMV),PSA 23203/0480

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)