

PathoScreen[®] Kit for AMV

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal with Anti-animal conjugate

Intended Use:

This ELISA is a qualitative serological assay for detection of AMV in plant foliage.

The test is proven to detect over 18 AMV isolates from all over the world.

Includes:

- Antibody-coated 96-well microtiter plates
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- RUB6 enzyme conjugate diluent, 1X ready to use
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- General extract buffer 3, GEB3
- [User Guide](#)



PathoScreen[®] Kit, Alfalfa mosaic virus (AMV), PSA 87601/0096

PathoScreen[®] Kit, Alfalfa mosaic virus (AMV), PSA 87601/0288

PathoScreen[®] Kit, Alfalfa mosaic virus (AMV), PSA 87601/0480

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com