

PathoScreen[®] Kit for PDV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Prune dwarf virus* (PDV) in stone fruit crops.

ELISA has proven to be an excellent diagnostic method for this pathogen as it reduces the time and labor involved with older techniques, such as grafting.

Includes:

- Antibody-coated 96-well microtiter plates
- Alkaline phosphatase enzyme conjugate
- RUB6 enzyme conjugate diluent, 1X ready to use
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- General extract buffer 3, GEB3
- [User Guide](#)



PathoScreen[®] Kit, Prune dwarf virus (PDV), PSA 98700/0096

PathoScreen[®] Kit, Prune dwarf virus (PDV), PSA 98700/0288

PathoScreen[®] Kit, Prune dwarf virus (PDV), PSA 98700/0480

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com