

## Bt-Cry1Ab/1Ac ELISA Kit

---

**Test Label:** Perox

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** NA

**Detection Reagent:** NA

### Intended Use:

---

This qualitative ELISA kit is used to detect the presence of Bt-Cry1Ab protein or the Bt-Cry1Ac protein expressed in transgenic crops. The test does not distinguish between Bt-Cry1Ab and Bt-Cry1Ac proteins.

This assay is suitable for testing both seed and leaf.

### Includes:

- Antibody-coated 96-well microtiter plates
- Peroxidase enzyme conjugate, concentrated
- RUB6 enzyme conjugate diluent
- TMB substrate solution
- Positive control
- PBST wash buffer, 20X concentrate
- [User Guide](#)



Bt-Cry1Ab/1Ac ELISA Kit, 96 testwells,PSP 06200/0096

Bt-Cry1Ab/1Ac ELISA Kit, 288 testwells,PSP 06200/0288

Bt-Cry1Ab/1Ac ELISA Kit, 480 testwells,PSP 06200/0480

Bt-Cry1Ab/1Ac ELISA Kit, 4800 testwells,PSP 06200/4800

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)