

## ELISA Reagent Set for Potato virus X (PVX)

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

### Intended Use:

This test can be used for detection of all strains of *Potato virus X (PVX)* in potato foliage, tuber sprouts and other solanaceous crops.

### Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for Potato virus X (PVX), 96 testwells,SRA 10000/0096

Reagent Set for Potato virus X (PVX), 500 testwells,SRA 10000/0500

Reagent Set for Potato virus X (PVX), 1000 testwells,SRA 10000/1000

Reagent Set for Potato virus X (PVX), 5000 testwells,SRA 10000/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)