

ELISA Reagent Set for TRSV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

Detection of all strains of *Tobacco ringspot virus* (TRSV) including Potato Calico strain. This test also detects *Potato black ringspot virus* (PBRV).

Includes:

- Capture antibody
- Alkaline phosphatase conjugated antibodies
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Tobacco ringspot virus (TRSV), 96 testwells,SRA 64001/0096

Reagent Set for Tobacco ringspot virus (TRSV), 500 testwells,SRA 64001/0500

Reagent Set for Tobacco ringspot virus (TRSV), 1000 testwells,SRA 64001/1000

Reagent Set for Tobacco ringspot virus (TRSV), 5000 testwells,SRA 64001/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com