

## ELISA Reagent Set for AMV

---

**Test Label:** Alkaline Phosphatase

**Test Format:** Triple Antibody Sandwich (TAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Monoclonal with Anti-animal conjugate

### Intended Use:

---

This ELISA is a qualitative serological assay for detection of *Alfalfa mosaic virus*(AMV) in plant foliage.

The test is proven to detect over 18 AMV isolates from all over the world.

### Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate
- RUB6 enzyme conjugate diluent, 1X ready to use
- General extract buffer 3, (GEB3)
- Tween 20
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for Alfalfa mosaic virus (AMV), 96 testwells,SRA 87601/0096

Reagent Set for Alfalfa mosaic virus (AMV), 500 testwells,SRA 87601/0500

Reagent Set for Alfalfa mosaic virus (AMV), 1000 testwells,SRA 87601/1000

Reagent Set for Alfalfa mosaic virus (AMV), 5000 testwells,SRA 87601/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)