

ELISA Reagent Set for Phyt

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

The *Phytophthora* ELISA test system is intended for use with plant samples exhibiting symptoms of *Phytophthora* infection. This test can also be used to test fungal culture samples.

The *Phytophthora* ELISA detects *Phytophthora* to the genus level and cannot differentiate species.

This test is the approved high throughput screening method in a U.S. national survey for *Phytophthora ramorum*, the causal agent of sudden oak death disease.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- RUB6, 1X Ready to use
- General extract buffer 2, GEB2 powder
- Postcoat 10
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Phytophthora (Phyt), 96 testwells, SRA 92601/0096

Reagent Set for Phytophthora (Phyt), 500 testwells, SRA 92601/0500

Reagent Set for Phytophthora (Phyt), 1000 testwells, SRA 92601/1000

Reagent Set for Phytophthora (Phyt), 5000 testwells, SRA 92601/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com