

Veratox for Hazelnut Allergen



Veratox for Hazelnut Allergen is used for the quantitative or qualitative analysis of hazelnut residues in food products such as cookies, crackers, chocolate bars, ice cream and cereals.


Features and Benefits

- Accurate and easy-to-use
- Quantitate or screen samples
- High specificity

Product Specifications

- Range of quantification: 2.5–25 ppm
- Controls provided: 0, 2.5, 5, 10 and 25 ppm hazelnut
- Testing time: 30 minutes
- Tests per kit: Up to 38

Results



Order Details

| Product # | Product Description |
|-----------|---|
| 8420 | Veratox for Hazelnut Allergen |
| 8429 | Allergen Extraction Kit - 20 samples |
| 8432S | Allergen Environmental Swabbing Kit – 100 samples |

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com