

# Veratox for Mustard Allergen



Veratox for Mustard Allergen is used for the quantitative analysis of mustard proteins in food products such as spices, sauces, dressings, meats and clean-in-place rinses.

## Features and Benefits

- Accurate and easy-to-use
- Quantitate or screen samples
- High specificity

## Product Specifications

- Range of quantification: 2.5–25 ppm
- Controls provided: 0, 2.5, 5, 10 and 25 ppm mustard
- Testing time: 30 minutes
- Tests per kit: Up to 38

## Results



## Order Details

Product #	Product Description
8400	Veratox for Mustard Allergen
8429	Allergen Extraction Kit - 20 samples
8432S	Allergen Environmental Swabbing kit – 100 samples

### India Contact:

## Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com), [www.atzlabs.com](http://www.atzlabs.com) ; [www.lifetechindia.com](http://www.lifetechindia.com)