

Whole Chromosome 4 Painting FISH Probe

Introduction

Whole Chromosome 4 Painting FISH Probe (WPP4) is designed to determine copy number and integrity of complete human chromosome 4 and detect abnormalities of the chromosome.

Intended Use

To measure copy number and integrity of human chromosome 4.

Cont.

WPP4 FISH Probe

Color

CytoOrange

Probe Design

Chr. 4



Whole Chromosome 4 Painting FISH Probe

Whole Chromosome 4 Painting FISH Probe, fragments with average length of ~300bp, covers entire human chromosome 4.

Not to Scale

Cat. No.

CT-WPP004-10-O

Volume

10 Tests (100 µL)

Signal Pattern Interpretation

Normal Patterns

20

Abnormal Patterns

Other Patterns

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com