

Product Information Sheet

B235

**Buffer Solution, pH Reference
Standard, pH 4.0 ± 0.02 at 25° C**

Properties

Form: Liquid
Appearance: Clear, Red
Application: pH Buffer Standard
Solubility:
Typical Working
Concentration:
Storage Temp: Room Temperature
Storage Temp of
Stock Solution:
Other Notes:

Application Notes

This is used for pH meter calibration

References

Revised 3/2007

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)
Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400
Email - customerservice@lifetechindia.com | customerservice@atzlabs.com