

AuPreP™ Nucleic Acid Purification Kits Range



AuPreP SPIN
Miniprep Kit
(Plasmid)



AuPreP Blood
Genomic DNA
Extraction Kit



AuPreP DNA easy
Plant Mini Kit



AuPreP PCR
Purification Kit



AuPreP Gel
Extraction Kit



AuPreP Plasmid
Maxi Kit



AuPreP Plant
RNA Maxi Kit



AuPreP RNA Easy
Midi Kit



AuPreP RNA
Extraction Mini
Kit



AuPreP RNA^m
Viral RNA
Extraction
Miniprep Kit