

# BioKits Walnut Assay Kit



BioKits Walnut Assay Kit is used for the detection and quantification of walnut in foods and environmental swabs.

## Product Specifications

- Lower limit of detection: 0.25 ppm walnut
- Range of quantitation: 2.4–120 ppm
- Testing time: 75 minutes
- Wells per kit: 48

## Results



## Order Details

Product #	Product Description
902085J	BioKits Walnut Assay Kit
901042J	BioKits Allergen Swabbing Kit – Contains 20 swabs, tubes and buffer

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com), [www.atzlabs.com](http://www.atzlabs.com) ; [www.lifetechindia.com](http://www.lifetechindia.com)