

ELISA Reagent Set for Prunus necrotic ringspot virus (PNRSV)

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Prunus necrotic ringspot virus* (PNRSV) in stone fruits.

Symptoms of PNRSV can include stunting, chlorosis, and leaf deformities. PNRSV is also the virus most associated with *Rose mosaic disease*.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Prunus necrotic ringspot virus (PNRSV), 96 testwells, SRA 30500/0096

Reagent Set for Prunus necrotic ringspot virus (PNRSV), 500 testwells, SRA 30500/0500

Reagent Set for Prunus necrotic ringspot virus (PNRSV), 1000 testwells, SRA 30500/1000

Reagent Set for Prunus necrotic ringspot virus (PNRSV), 5000 testwells, SRA 30500/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com