

ELISA Reagent Set for Beet western yellows virus (BWYV)

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

This ELISA test is a qualitative serological assay for detection of *Beet western yellows virus* (BWYV) in vegetable, floriculture, and legume crops.

This test has shown cross-reaction to BYDV-rpv isolates.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Beet western yellows virus (BWYV), 96 testwells, SRA 32100/0096

Reagent Set for Beet western yellows virus (BWYV), 500 testwells, SRA 32100/0500

Reagent Set for Beet western yellows virus (BWYV), 1000 testwells, SRA 32100/1000

Reagent Set for Beet western yellows virus (BWYV), 5000 testwells, SRA 32100/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com