

ELISA Reagent Set for Tomato spotted wilt virus (TSWV)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Monoclonal

Detection Reagent: Monoclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Tomato spotted wilt virus* (TSWV) in plant foliage.

The test is proven to detect over 38 TSWV isolates from all over the world.

This test has been found to exhibit cross-reactivity with some isolates of *Groundnut ringspot virus* (GRSV) and *Tomato chlorotic spot virus* (TCSV).

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Tomato spotted wilt virus (TSWV), 96 testwells,SRA 39300/0096

Reagent Set for Tomato spotted wilt virus (TSWV), 500 testwells,SRA 39300/0500

Reagent Set for Tomato spotted wilt virus (TSWV), 1000 testwells,SRA 39300/1000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com



Life Technologies™

Your Molecular & Cell Technology Partner



agdia®

Reagent Set for Tomato spotted wilt virus (TSWV), 5000 testwells, SRA 39300/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com