

## ELISA Reagent Set for TBSV

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal


### Intended Use:

This ELISA is a qualitative serological assay for detection of TBSV in plant foliage.

This test is also known to cross-react with *Lettuce necrotic stunt virus* (LNSV), *Moroccan pepper virus* (MPV) and *Pelargonium leaf curl virus* (PLCV)

### Includes:

- Capture antibody
- Alkaline phosphatase conjugated antibodies
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for Tomato bushy stunt virus (TBSV), 96 testwells, SRA 45400/0096

Reagent Set for Tomato bushy stunt virus (TBSV), 500 testwells, SRA 45400/0500

Reagent Set for Tomato bushy stunt virus (TBSV), 1000 testwells, SRA 45400/1000

Reagent Set for Tomato bushy stunt virus (TBSV), 5000 testwells, SRA 45400/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)