

ELISA Reagent Set for Pstew

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal with Anti-animal conjugate

Intended Use:

This ELISA is a qualitative serological assay intended for detection of *Pantoea stewartii* bacteria the causal agent of Stewart's wilt.

Agdia's ELISA test for *Pantoea stewartii* provides an effective, lab-based method for phytosanitary assessment.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for *Pantoea stewartii* (Pstew), 96 testwells, SRA 52000/0096

Reagent Set for *Pantoea stewartii* (Pstew), 500 testwells, SRA 52000/0500

Reagent Set for *Pantoea stewartii* (Pstew), 1000 testwells, SRA 52000/1000

Reagent Set for *Pantoea stewartii* (Pstew), 5000 testwells, SRA 52000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com