

ELISA Reagent Set for IYSV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Iris yellow spot virus*(IYSV) in plant foliage.

The test is proven to detect over 6 IYSV isolates from various regions throughout the world.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Regent Set, Iris yellow spot virus (IYSV), 96 testwells,SRA 60500/0096

Regent Set, Iris yellow spot virus (IYSV), 500 testwells,SRA 60500/0500

Regent Set, Iris yellow spot virus (IYSV), 1000 testwells,SRA 60500/1000

Regent Set, Iris yellow spot virus (IYSV), 5000 testwells,SRA 60500/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com