

ELISA Reagent Set for KLV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This ELISA test is a qualitative serological assay for detection of *Kalanchoe latent virus* (KLV) in Kalanchoes. Two strains of the virus, KaV-1 and KaV-2, have been described.

Since KLV often does not show visible symptoms, infection is difficult to determine without a diagnostic laboratory technique.

Includes:

- Capture antibody
- Alkaline phosphatase conjugated antibodies
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Kalanchoe latent virus (KLV), 96 testwells, SRA 62000/0096

Reagent Set for Kalanchoe latent virus (KLV), 500 testwells, SRA 62000/0500

Reagent Set for Kalanchoe latent virus (KLV), 1000 testwells, SRA 62000/1000

Reagent Set for Kalanchoe latent virus (KLV), 5000 testwells, SRA 62000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com