

ELISA Reagent Set for ZYMV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This is a qualitative serological assay for detection of *Zucchini yellow mosaic virus*(ZYMV) in plant foliage.

The test has been used for crop screens by Agdia Testing Services Department using cucurbit hosts, including squash, cucumber and melon.

The ZYMV ELISA detects the most important isolates from all over the world, including CT, USDA, SJBCA, CA, IT, NY, FL, Z18

Includes:

- Capture antibody
- Alkaline phosphatase conjugated antibodies
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Zucchini yellow mosaic virus (ZYMV), 96 testwells, SRA 77700/0096

Reagent Set for Zucchini yellow mosaic virus (ZYMV), 500 testwells, SRA 77700/0500

Reagent Set for Zucchini yellow mosaic virus (ZYMV), 1000 testwells, SRA 77700/1000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com



Reagent Set for Zucchini yellow mosaic virus (ZYMV), 5000 testwells,SRA 77700/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com