

ELISA Reagent Set for PRMV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Peach rosette mosaic virus* (PRMV) in peach and grape crops.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Peach rosette mosaic virus (PRMV), 96 testwells, SRA 87700/0096

Reagent Set for Peach rosette mosaic virus (PRMV), 500 testwells, SRA 87700/0500

Reagent Set for Peach rosette mosaic virus (PRMV), 1000 testwells, SRA 87700/1000

Reagent Set for Peach rosette mosaic virus (PRMV), 5000 testwells, SRA 87700/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com