

Water DNA Isolation

UltraClean™ Water DNA Isolation Kit (0.45 um)



Other Products


- UltraClean™ Water DNA Isolation Kit (0.45 um)
- UltraClean™ Water DNA Isolation Kit (0.22 um)
- UltraClean™ Water DNA Isolation Kit (No filters)
- Water Filters
- UltraClean™ Water DNA Isolation Kit Components

For isolation of DNA from microorganisms in small or large volume water samples.


Product Features and Benefits

- ☐ **Rapid Protocol** - Purify DNA from microorganisms in water in 35 minutes.
- ☐ **High Quality DNA** - Yields DNA greater than 50 kb in length ready for PCR, restriction digests and other downstream applications.
- ☐ **Flexible Format** - Useful for small or large volumes of water.
- ☐ **Complete** - All tubes and reagents required are included in kit and ready to use.
- ☐ **Safe** - No hazardous organic solvents or enzymes required.

DNA isolated from sewage treatment plant in Carlsbad, California using the UltraClean™ Water DNA Isolation Kit



Lane 1: 1 kb Marker
Lanes 2, 3: Primary treated sewage
Lanes 4, 5: Secondary treated sewage



To order contact:

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400 Fax: +91-11-42208444

Email: customerservice@atzlabs.com

Web: www.atzlabs.com