

Product Information Sheet

O139

Orchid Maintenance/Replate Medium Without Charcoal and Agar

Properties

Form:	Powder
Appearance:	White to Yellow powder
Application:	Orchid Culture
Solubility:	Water
Typical Working Concentration:	25.31 g/L
Storage Temp:	2 – 6° C
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes:	pH = 4.75 – 5.75

Formula

(mg/L)

Ammonium Nitrate	825
Boric Acid	3.1
Calcium Chloride, Anhydrous	166
Cobalt Chloride•6H ₂ O	0.0125
Cupric Sulfate•5H ₂ O	0.0125
Na ₂ EDTA•2H ₂ O	37.3
Ferrous Sulfate•7H ₂ O	27.85
Magnesium Sulfate, Anhydrous	90.35
Manganese Sulfate•H ₂ O	8.45
Molybdic Acid (Sodium Salt) •2H ₂ O	0.125
Potassium Iodide	0.415

Potassium Nitrate	950
Potassium Phosphate, Monobasic	85
Zinc Sulfate•7H ₂ O	5.3
MES (Free Acid)	1000
myo-Inositol	100
Nicotinic Acid (Free Acid)	1
Peptone from Meat	2000
Pyridoxine•HCl	1
Sucrose	20,000
Thiamine•HCl	10

Application Notes

Plant Tissue Culture Tested

Plant species: Many epiphytic orchid species

This medium is a modification of our original orchid replate medium (P668) except without charcoal (and agar). This medium permits the culture of orchids in a liquid system or for orchid species that do not require charcoal to neutralize phenolics produced by the plants.

Revised 12/2012

India Contact

Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, Road 44, Pitampura, Delhi - 110034 (India)

Tel: +91-11-4220-8000; 4220-8111; 4220-8222 Fax: +91-11-4220-8444, Mobile: +91-98105-21400

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com